
elexxion

powered by technology

An abstract graphic featuring flowing, wavy lines in shades of blue and white, set against a dark background. The lines create a sense of motion and depth, resembling a liquid surface or a digital wave.

miscellaneous application studies

Table of contents

Title	Author	Date of Expiry	Page
Depigmentation of gingivae and lip with digital pulsed diode laser – an integral part of cosmetic dentistry	Dr. Kenneth Luk, Hong Kong	2009	5
Der Diodenlaser als Hard- und Therapielaser	Frank W. Förster	2008	9
Der Laser – ein rückgekoppelter optischer Verstärker in Selbsterregung	Prof. Dr. Axel Donges	2006	17
Clinical application of a digital pulsed diode laser in depigmentation therapy	Dr. Kenneth Luk	2005	25
Ein Jahr Elexxion claros	Dr. med. dent. Michel Vock		31
The evolution of the diode laser	Olaf Schäfer		35

Depigmentation of gingivae and lip with **Digital Pulsed Diode Laser**—an integral part of Cosmetic Dentistry

author_Kenneth Luk, Hong Kong

Melanin is the most common natural pigment contributing to normal colour of the gums. Although ethnic and physiologic melanin pigmentation is not a medical pathosis, the appearance of pigmented gums is considered unaesthetic and overlooked at times in cosmetic dentistry assessments. Various laser wavelengths have been reported for removal of melanin pigment with good results. Removal of mucosa to the basal layer is the common procedure using lasers or conventional abrasive technique using diamond burs. Pigmentation of the lip is another aesthetic concern to patients affecting their confidence and quality of life.

Background

Diode laser at 810 nm is transmissive through water, and is attracted to pigment such as haemoglobin and melanin. The continuous wave or the mechanically gated pulsed mode both produce a long tissue interaction time, which require much longer relaxation time for tissue to cool down. Hence, low power (1 to 3 W) is recommended for most soft tissue procedures. The use of "High Fluency Technique", for 980 nm diode wavelength developed by Dr Mick Swick, is able to increase power and reduce collateral tissue damage and heat stacking effect by using water irrigation as a cool-

Case 1

ant. The major precaution for any soft tissue laser procedure is to avoid deep penetration of the laser energy, which can cause collateral thermal damage to the underlying bone or adjacent tooth structures, especially when using a non-contact and/or non-initiated tip. Digital pulsed diode laser is the 4th generation of diode laser technology. Dr Claus Neckel compared fibromas excised by the pulsed diode laser with continuous wave diode laser. He demonstrated the substantial reduction in collateral thermal damage histologically with high power (30 W), short pulse (10 μs) with high frequency (10,000 Hz). As melanin is deposited at the basal layer of the

Fig. 1 Pre-op.
Fig. 2 Immediate post-op.
Fig. 3 Pre-op.
Fig. 4 18 months post-op.
Fig. 5 4.5 year post-op.

case_report

mucosa, removal of this pigment by various lasers have been reported by direct ablation (cutting) of the mucosa to this layer. Energy delivery time ranges from 10 to 30 minutes in an area of 1st premolar to 1st premolar of one arch. The technique described in this case series differs from other laser techniques in terms of tissue absorption and procedure time.

Procedure for depigmentation of gingivae

Elxion Claros DPL Laser set at 30 W, 16 μ s, 20,000 Hz; average power close to 10 W.

Under local anaesthesia, a non-initiated 600 μ m fiber delivers the energy at a distance of 2 mm to 5 mm to the pigmented area with constant movement under water irrigation as coolant. There is no direct ablation of the pigmented mucosa but rather the hemoglobin and melanin absorbing the laser energy. Procedure time ranges from two minutes to seven minutes in an area of 1st premolar to 1st premolar of one arch.

Case 1

The procedure took approximately 2 minutes to complete between upper left and right premolar region. 4.5 year post-op showed mild relapse of pigmentation but patient is still happy with the colour. Laser soft tissue crown lengthening was carried out on upper anterior incisors and canines.

Case 2

The procedure took approximately 7 minutes to complete between upper left and right premolar region.

Case 2

Fig. 1_Pre-op.
Fig. 2_1.5 year post-op.

Case 3

Fig. 1_Pre-op.
Fig. 2_3 months post-op.


Table 1_Graphic display of chromophores of various wavelengths. The 810 nm diode is indicated.

Case 3

The procedure took approximately 2.5 minutes to complete between lower left and right premolar region. Ceramic crowns replacing all porcelain fused to metal crowns

Procedure for depigmentation of lip

The procedure is the same as for gingival depigmentation. As it took only seconds to complete the procedure, there was no need for water cooling.

Case 4

The procedure took 8 seconds to complete under local anaesthesia.

Case 5

Procedure time: 3 separate visits under topical anaesthetic over 8 months (4 sec, 5 sec, 11sec)

_Case 4


Fig. 1_Pre-op.
Fig. 2_Immediate post-op.
Fig. 3_1.5 year post-op.

_Case 5

Fig. 1_Pre-op.
Fig. 2_6 months post-op.

_Results

There are two possible immediate visual post-op results. In ideal case (Case 1), the mucosa turns pink without any signs of surface mucosal damage. Sub surface coagulation of blood vessels gave a pink coloured appearance. No laser peeling of mucosa was noted post-operative result. Other cases vary in degree of immediate aesthetic improvement or darkening of the pigmented mucosa where absorption by melanin is the dominant. Laser peeling of the mucosa were reported within a few days of treatment. In all cases, there were zero to mild post-op pain. No analgesics were required.

_Conclusion


This technique requires specific high power settings. Operator should be aware of tissue interaction during the procedure and adapt accordingly. Long term relapse of the condition is comparable to direct removal of pigment by lasers or other techniques. Pa-

tients were all satisfied with the results. The main advantages with this technique are minimal invasive technique, short procedure time and immediate aesthetic results in some cases.

Advantages	Disadvantages
Minimal invasive technique	Local anaesthetics
Short treatment time	Immediate aesthetic results cannot be achieved in every case
Pain free post-op (Immediate aesthetic results)	Technique sensitive
	Specific power settings not applicable to other diode lasers

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Kenneth Luk, BDS, DGD
Hong Kong

Kenneth Luk, BDS,
DGD

CLINICAL APPLICATION OF A DIGITAL PULSED DIODE LASER IN DEPIGMENTATION THERAPY

Dr. Kenneth Luk received his BDS degree from the University of Liverpool and was subsequently awarded the Diploma in General Dental Practice from The Royal College of Surgeons (England). Currently, he is a part-time lecturer at The University of Hong Kong and also maintains a private practice in Hong Kong. He has achieved fellowship status from the International Congress of Oral Implantologists (ICOI). Dr. Luk has incorporated the use of lasers in his practice since 2003. He is a member of the Academy of Laser Dentistry (ALD) and has achieved Standard Proficiency in various laser wavelengths. He also serves on ALD's International Relations Committee.

Disclosure: Dr. Luk is a presenter for Elexion and receives an honorarium for his services.

Pretreatment

A. Case Outline

A 26-year-old female of Chinese ancestry presented with a chief concern of unaesthetic dark-colored gums. She noted this as a cosmetic problem, present for more than 10 years. She was healthy with a negative medical history and smokes 4 cigarettes per day.

Her skin complexion is fair (Figure 1). She was not a mouth breather and her lips were competent. Dark, unaesthetic gingival tissue was noticeable in her smile profile (Figure 2).

Intraoral examination revealed melanin pigmentation on both upper and lower labial gingival tissues (Figure 3). There were no restorations in her dentition and the periodontal condition was good.

B. Diagnosis

1. Provisional Diagnosis

Congenital melanin pigmentation of the labial gingivae.

2. Treatment Plan

The proposed plan was to use laser energy to conservatively ablate the pigmented gingival epithelium. Two laser wavelengths were to be employed in sequence, with the use of the second wavelength dependent on the outcome of the first.

The patient was informed of possible surface tissue ablation on the treated sites, and was told that relapse of the condition is possible.

The patient agreed to the plan after an informed consent/refusal discussion took place.

The lower anterior soft tissue (adjacent to teeth #22-27) was treated once with an erbium:YAG laser (Opus Duo, OpusDent, Yokneam, Israel), with an emission wavelength of 2940 nm. The operating parameters were 400 mJ, 10 Hz, and a 1000-micron sapphire tip was used. A protective gingival spacer was also used (Figure 4). The procedure was performed under topical anaesthesia with no water spray until the target pigment area was removed. Total exposure duration was estimated at 1 minute and 40 seconds. Figures 5 and 6 show the immediately postoperative result. The hemostasis was very good but the patient reported some discomfort during the procedure. Six weeks later, the patient returned for a follow-up examination and there was no significant reduction in the pigmentation (Figure 7).

The patient then asked whether the diode laser would be effective for depigmentation, since her outcome experience was not good with the erbium laser. She was then offered to have the procedure repeated, this time using a diode laser on the other pigmented area.

Depigmentation of upper gingivae (adjacent to teeth #5-12) under local anesthetic with a digital pulsed diode laser (810 nm wavelength) was then agreed upon, and is the focus of this manuscript.

3. Indications, Considerations, and Precautions for Lasers

Various laser wavelengths have shown success in melanin depigmentation of gingiva (see page 16).

The graph (Diagram 1) shows the relative absorption of different laser wavelengths. The erbium:YAG (2940 nm) and CO₂ (10,640 nm) lasers produce soft

tissue surface ablation that is relatively non-selective because of the high water absorption as well as minimal tissue penetration. The dental erbium laser, with its microsecond pulse width and low repetition rate, typically produces less hemostasis than the CO₂ laser.

The Nd:YAG laser at 1064 nm is a versatile soft tissue laser. Similar to diode lasers, it reacts with the chromophores of melanin and hemoglobin. The Nd:YAG, transmissive through water, can penetrate into soft tissue. Its free-running pulse mode allows for thermal tissue cooling between the short pulses.

The diode lasers at 810 and 980 nm are also transmissive through water, and are attracted to pigmented tissue. The continuous wave or the mechanically gated pulsed mode both produce a long tissue interaction time. This requires much longer relaxation time for tissue to cool (Diagram 2). Hence, low power (1 to 3 W) is recommended for most soft tissue procedures. The use of a "high fluency technique" for the 980 nm diode wavelength, developed by Dr. Mike Swick, enables the clinician to increase power and reduce potential collateral tissue damage from the heat-stacking effect by using water irrigation as a coolant.

A major precaution during any soft tissue laser procedure is to avoid collateral thermal damage to underlying bone or adjacent tooth structures, especially when using a noncontact and/or noninitiated tip.

The digital pulsed diode laser is the 4th generation of diode laser technology. The aim of using a digital pulsed diode is to achieve high power intensity and short pulse width (16 µsec) with penetration only into the pigmented area, thereby minimizing the surface ablation effect and

the heat-stacking effect beyond the target tissue (Diagram 3).

Treatment

A. Laser Operating Parameters

Digital Pulsed Diode Laser (Claros, Elexxon, Radolfzell, Germany)

Wavelength: 810 nm

Delivery System: Fiberoptic

Beam: 600 μ m noninitiated

Power: 30 W

Tissue Exposure Mode: Digital Pulsed 16 μ sec

Pulse energy: 0.48 mJ

Pulse Rate: 20,000 Hz

Average power: 9.99 W

Exposure duration: 2 minutes

B. Laser Safety Check

A laser safety check was carried out.

The operating room was closed with a laser-in-use sign posted outside.

All personnel and the patient were issued protective glasses (L6). This is equivalent to an optical density of 6. Note: Normal protective glasses (L4) from other laser devices for the same wavelength are not recommended for this laser device.

As an added safety measure, a remote controlled signal light (battery-operated flashlight and a transmitter) is placed outside the operating room. When the laser is firing, the signal light flashes with an orange-colored light every second.

C. Treatment Delivery Sequence

Buccal local infiltration was given using 2 carpules of 2% lignocaine with 1:80000 adrenaline (Xylocaine, Astra, North Ryde, Australia).

OpalDam® (Ultradent Products, South Jordan, Utah) was placed on the cervical margin from tooth #5 to 12 to protect the teeth from laser exposure.

The 600-micron noninitiated fiber was placed 2 mm from the target tissue of teeth #5, 6, and 7. The fiber was directed at right angles to the target tissue. The fiber was kept in motion at all times, and the total exposure time was 26 seconds. Tissue blanching and instant depigmentation were visible. No surface ablation was noted.

The decision was then made to add water irrigation to the procedure for the

remaining soft tissue areas (teeth #8-12). The aim was to avoid dehydration by high-speed evacuation by keeping the surface tissue cool and moist. The same parameters as above were used for an additional exposure time of 1 min 14 sec.

As noted, the total laser exposure lasted 1 minute and 40 seconds, and as a result of the water irrigation, there was no sign of tissue blanching.

D. Postoperative Instructions

Postoperative instructions were given to the patient as follows:

- Diet control: Avoid colored food additives like soya sauce, and avoid spicy food for one week.
- Refrain from smoking for one week.
- The patient was given a prescription for ibuprofen if needed for discomfort.
- Rinse with warm salt water and maintain good oral hygiene.

E. Assessment of Treatment Outcome

There was no surface tissue damage and instant depigmentation was visible.

There were no intraoperative or postoperative complications, and the prognosis was very good.

Follow-Up Care

A. Assessment of Treatment

When the patient returned to her office immediately after treatment, her colleagues noticed the marked difference in her smile. She reported mild discomfort after the anesthetic wore off, although she took no analgesics. She also reported that her teeth were sensitive at times on the first day.

The patient returned for her first follow-up 6 days following treatment (Figure 8). The tissue was pink in color, and an overexposed area around tooth #7 was healing within normal limits.

At 3 weeks, the patient returned for inspection and there was no darkening in color of the gingivae (Figure 9).

At the 7-week follow-up visit, the gingival tissue had not darkened and in fact remained lighter in color than the preoperative condition. The teeth were free of symptoms (Figure 10).

At 6 months, the patient noticed mild

darkening of the gingivae. She was still very satisfied with the result (Figure 11). Electric pulp testing (EPT) and cold tests were within normal limits. Periapical radiographs revealed no sign of any apical lesions or abnormalities (Figures 12-13).

B. Conclusion

The result of this technique is comparable to other techniques used by other laser wavelengths. The immediate postoperative cosmetic appearance of the gingivae is an advantage with this procedure which the author named "Atraumatic Depigmentation Therapy."

Since the treatment was performed, the patient has been more confident with her smile (Figures 14-16). Although local anesthetic was required for treatment with the diode laser, this technique was preferred to the other procedure performed on the lower gingivae, both in terms of less discomfort during the procedure and improved gingival appearance immediately postoperative. The patient's lower gingivae did not show when smiling, and she was not concerned with further treatment of that tissue. She was very satisfied and expressed a desire to whiten her teeth.

The digital pulsed diode laser enhances previous diode technologies which have a longer gated pulse width. It thus offers more energy delivery options for desired tissue interactions. Although this case gives successful cosmetic results, more cases and longer postoperative follow-ups should be carried out to give a comprehensive assessment of this technique.


Diagram 1: Graphic display of chromophores of various laser wavelengths. The 810 nm diode laser is indicated with a dotted line.


Diagram 2: Graphic display of the constant tissue temperature rise when exposed to continuous wave laser energy.


Diagram 3: Graphic display of tissue temperature rise and fall in response to the small interaction time of a short-pulsed laser.

Figure 1: Preoperative, frontal view

Figure 2: Preoperative smile

Figure 3: Preoperative anterior gingival tissue showing hyperpigmentation

Figure 4: Gingival spacer shield fitted to the Er:YAG laser handpiece and tip

Figure 5: Immediate postoperative view of the mandibular right side, Er:YAG laser used

Figure 6: Immediate postoperative view of the mandibular left side, Er:YAG laser used

Figure 7: Six weeks postoperative, Er:YAG laser

Figure 8: Six days postoperative, maxillary anterior, diode laser used

Figure 9: Three weeks postoperative, maxillary anterior, diode laser used

Figure 10: Seven weeks postoperative, maxillary anterior, diode laser used

Figure 11: Six months postoperative, maxillary anterior, diode laser used

Figure 12: Six-month postoperative radiograph, maxillary anterior right side

Figure 13: Six-month postoperative radiograph, maxillary anterior left side

Figure 14: Six months postoperative, right profile smile

Figure 15: Six months postoperative, left profile smile

Figure 16: Six months postoperative, frontal view

Lasers and Depigmentation

Different laser wavelengths have been used for removal and reduction of congenital gingival pigmentation. Other techniques include conventional periodontal gingivectomy and gingivoplasty procedures, as well as the use of cryogenics. A sampling of the literature reveals several laser references with successful outcomes:

1. Yousuf A, Hossain M, Nakamura Y, Yamada Y, Kinoshita J, Matsumoto K. Removal of gingival melanin pigmentation with the semiconductor diode laser: A case report. *J Clin Laser Med Surg* 2000;18(5):263-266.
2. Atsawasuwan P, Greethong K, Nimmanon V. Treatment of gingival hyperpigmentation for esthetic purposes by Nd:YAG laser: Report of 4 cases. *J Periodontol* 2000;71(2):315-321.
3. Tal H, Oegiesser D, Tal M. Gingival depigmentation by erbium:YAG laser: Clinical observations and patient responses. *J Periodontol* 2003; 74(11):1660-1667.
4. Esen E, Haytac MC, Oz IA, Erdogan O, Karsli ED. Gingival melanin pigmentation and its treatment with the CO₂ laser. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2004;98(5):522-527.

In the previous issue of the *Journal of the Academy of Laser Dentistry*, Dr. Kenichi Takahashi described his treatment of hyperpigmentation using an Er:YAG laser in a contact mode with an 800-micron diameter fiber and an average power of 2 W (200 mJ, 10 Hz), which is a different parameter and spot size from Dr. Luk's technique. Dr. Takahashi had good results, although the reader will note that surface ablation did occur and hemostasis was only adequately managed with the erbium wavelength.

In *Wavelengths*, Dr. Michael Swick explains his pioneering technique of using a copious water spray with high fluence from the 980 nm diode laser. His clinical results show one good solution to the problem of excessive thermal energy effects.

1. Swick MD. Cosmetic diode laser gingivectomy with frenectomy. *Wavelengths* 2000;8(4):19.
2. Swick MD. A char-free technique for the Ceralas D15 diode laser. *Wavelengths* 2000;8(4):20.

Dr. Luk's instrumentation demonstrates another method of attenuating the tissue temperature rise by using very short pulses of 810 nm diode laser energy.

As always, the laser energy must be carefully directed toward the target tissue to remove the pigmented gingiva; each wavelength will have unique interaction properties, and the clinician must carefully evaluate the progress of the surgery.

Lasers in Dentistry

In this issue, we feature two laser wavelengths for a variety of soft tissue and hard tissue procedures

- Dr Kenneth Luk uses a diode laser with digitally pulsed technology to remove melanin hyperpigmentation from otherwise healthy gingival tissue. His case shows a comparison between his use of an erbium laser and then his subsequent use of a diode laser.

The *Journal* is again pleased to publish clinical cases from three of the recent successful Advanced Proficiency candidates who presented during the ALD's 12th Annual Conference in New Orleans.

- Dr. Emil Litvak performs removal and recontouring of hyperplastic gingival tissue that occurred following orthodontic treatment. He carefully employs an Er:YAG laser to recreate a more normal amount of periodontal soft tissue.
- Dr. Glenda Payas has two clinical case studies using an Er:YAG laser. In the first, she performs soft tissue crown lengthening on the upper anterior arch for a patient with hyperplastic gingival tissue; the second case shows removal of several bony protuberances that occurred following full-mouth extractions. She employs this wavelength to its full advantage in safely treating gingiva and bone.
- Dr. Ronald Schalter also presents two clinical case studies with an Er,Cr:YSGG laser. In the first case, he prepares the upper and lower anterior teeth for direct composite restorations. His second case demonstrates a fibroma excision. In both instances, Dr Schalter demonstrates the usefulness of this wavelength for hard and soft tissue treatments.

These cases show the numerous clinical indications for use of both the diode and the erbium family of lasers and once again demonstrate how dental lasers continue to be utilized for good patient care.

Case Studies

- **Clinical Application of a Digital Pulsed Diode Laser in Depigmentation Therapy**
KENNETH LUK, BDS, DGDP – HONG KONG
- **Er:YAG Laser-Assisted Gingivectomy after Orthodontic Treatment**
EMIL LITVAK, DDS – YAHUD, ISRAEL
- **Clinical Application of an Er:YAG Laser in Gingival Recontouring of Maxillary Teeth**
- **Clinical Applications of an Er:YAG Laser in Removal of Bony Protuberances Following Extractions and Mandibular Tori Removal**
GLENDA PAYAS, DMD – TULSA, OKLAHOMA
- **Multiple Restorative Treatments Using an Er,Cr:YSGG Laser**
- **Fibroma Removal Using an Er,Cr:YSGG Laser**
RONALD W. SCHALTER, DDS – HUDSON, MICHIGAN

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Literatur

DÖRTBUDAK O (2001): Photodynamic therapy for bacterial reduction of periodontal microorganisms. Journal of oral Laser Applications 1: 115-118

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MORITZ A. et al. : Bacterial reduction in periodontal pockets through irradiation with a diode laser. Journal of clinical laser medicine 1/97, 33-37 (1997)

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

The evolution of the diode laser

In modern dentistry applications, two different types of laser occupy dominant positions: gas lasers and solid-state lasers. The first of these categories includes the Co₂ laser (10 600 nm) and the argon laser (480 + 515 nm), while examples of the solid-state laser are the Nd:YAG (1064 nm), Er:YAG (2940 nm), Er:CrYSGG (2780 nm) and the diode laser (mainly 810 and 980 nm). As the name suggests, gas lasers use gas mixtures as their active laser medium whereas solid-state lasers usually employ a crystal as their active laser medium. In classic solid-state lasers, the laser medium takes the form of an optical crystal which is excited by a flashlight (a process known as “pumping”). The laser crystal and flashlight are located in a unit known as a resonator, a chamber comprising 2 mirrors between which the laser process takes place. At a low percentage rate, one of these mirrors is permeable to laser light, enabling usable light energy to be disengaged from the process. A fundamental property of this design principle involves a certain level of laser energy being stored in the resonator which can be released in a very short pulse. The shorter the energy release time (i.e. the pulse time), the lower the thermal impact of the laser beam. In eye surgery applications, tissue can be removed by very short pulses of laser light to correct vision defects, and this can be done without damaging this tissue with heat. In the world of dentistry, manufacturers of Er:YAG lasers are committed to achieving the shortest possible pulse times. This increases the removal rate for hard tooth material while at the same time reducing the sensation of pain.

Bildertext: *Electrical energy* *Flashlight*
Energy supply, pumping
Mirror *Resonator* *Mirror*

Semiconductors or diode lasers occupy a special position among solid-state lasers. These do not have a resonator in any conventional sense of the term because the laser medium takes the form of a semiconductor chip. In keeping with all semiconductors, the diode laser has a p-layer and an n-layer. The laser beam is generated precisely in the boundary layer between the p-layer and n-layer. The great advantage of this is that you do not need any flashlights for optical pumping. Instead, the laser element can be excited directly with electrical current. This gives rise to a very high efficiency rating of approx. 35 % (comparison with YAG laser: approx. 1-5 %). For this reason, diode lasers do not require large cooling units and can be built in relatively compact, lightweight units.

Bildertext

Fig. 15 Diagram of a laser diode

A minor disadvantage of the diode laser is that no meaningful levels of energy can be stored due to the absence of a resonator chamber. This means that the laser can only be “on” or “off”. For this reason, diode lasers are typically operated in continuous or CW (continuous wave) mode.

Diode laser technology first started being used in the early 1980’s in the context of what became known as therapy lasers. They were also known as soft lasers or LLLT lasers (low level laser therapy). The spectrum ranges from 1 mW at 635 nm – equivalent to a standard laser pointer – up to 200 mW at 810 nm. Everyone is free to form their own judgement about the possible differences in operating method.

Around about 1995, the first diode laser devices made their appearance, with a power rating of roughly 6 W at 810 nm. For the first time, it was now possible to utilise the thermal effect of diode laser radiation. Numerous clinical studies proved that superlative results could be achieved with diode laser radiation in the 1.0 Watt range when applied to the decontamination of periodontal abscesses and to root canal work. This is due to the good absorption characteristics of the laser light in gram-negative and positive bacterial cells at this wavelength. The accompanying high absorption rate for haemoglobin observed means that, at higher energy levels, surgical incisions can be made using these systems, accompanied by an almost complete absence of bleeding. As time went by, other meaningful possible applications were added and these are now considered to be clinically secure. As a typical example, we would point to the treatment for periodontal implants (“periimplantitis”) which has been superbly documented by a wide range of universities around the world.

To further improve performance in the surgical sector, it was appropriate to boost the output power of the diode laser. Initial tests with 10 and 15 Watt diodes demonstrated as early as 1999 that success could not be achieved simply by increasing the light energy. A power rating of about 3-4 Watts CW led fairly rapidly to carbonisation of the tissue. Due to the resultant toxic by-products, this effect is therefore undesirable. In addition, wound healing and the sensation of pain are both adversely affected. The cause of this carbonisation is not therefore related to wavelength, instead it can be traced to the timing characteristics of the laser. Since the laser is continuously enabled, this gives rise to a continuous increase in tissue temperature. The burned


tissue also exhibits a higher absorption coefficient than healthy tissue. As a consequence, this rise in temperature is disproportionately high.

If you wish to further increase the cutting speed – without damaging the tissue through severe carbonisation, there is only one way left open: simulation of the pulse characteristics of a solid-state laser, e.g. the Er:YAG. Since the diode laser source, as described at the outset, is not able to “store” energy in any significant amounts, the laser diode needs to be able at all times to generate the desired pulse power. The CW performance of the laser diode therefore corresponds to maximum pulse power.

Initial tests were conducted in 2000 with an 810 nm diode laser, an output rating of 20 W and minimum pulse times of approximately 50 microseconds (millionths of a second). In this process, the time between 2 pulses was set to between 100 and 450 microseconds, equivalent to a mean power rating of 2.0 – 6.7 W.

The histological and practical results were very promising indeed. Compared with a CW diode laser of the same power rating, the necrosis zone in the cut section was almost halved while, at the same time, cutting speed was increased substantially.

Based on these results, work commenced in 2002 on the development of a 30 W diode laser (Elexxion) with a minimum pulse duration of just 9 microseconds (millionths of a second). Up to 20 000 pulses per second can be issued across a variable range and, at appropriate pulse widths, the calculated mean power rating was restricted to a maximum of 10 watts. Technical implementation of these short pulses at the high pulse power rating of 30 W necessitated the development of a completely new digital electronic activation system. This type of pulsing process is also referred to as the digital pulsed laser (DPL[®]) technique. Furthermore, fibre applicators with a special coating had to be developed, capable of withstanding these loadings. Nevertheless, it is still virtually impossible to cause tissue burns, provided that the power settings defined by Elexxion (the manufacturer) are complied with.

Bildertext

Standard pulse technology
CW diode

Clinical trials have demonstrated that a cutting speed in soft tissue can be achieved which actually outperforms a Co₂ laser. At the same time, this laser light cuts more gently than ever before. Previous objections to diode laser surgery (“cutting with a hot glass needle”) can now be refuted very credibly indeed.

Something which works so effectively in surgical situations ought to be equally effective in paradontology and endodontics. For the last ten years, these sectors have been using diode lasers with power ratings of 1.0 W for decontamination work. This restricted power rating is chosen to prevent undesirable thermal effects from arising. It is conceivable, through the use of extremely short 30 W pulses, to improve effectiveness substantially without entailed undesirable side effects. Initial clinical trials are now in progress.

We shall see what the future holds. In 5 years time, we might already be working with 50 or even 100 W diode lasers. One thing is for certain though: the future of the diode laser is already in its infancy.

Olaf Schäfer
Chief Executive
elxxion GmbH Radolfzell

elexxion AG

Schützenstraße 84 · 78315 Radolfzell · Germany
Tel. 0049 (0) 7732-822 99 0 · Fax 0049 (0) 7732-822 99 77
info@elexxion.com · www.elexxion.com