

Comparison of Diode Lasers in soft- tissue surgery using
cw- and superpulsed mode: an in vivo study.

by

Dr. Ralf Borchers

for the degree of
Master of Science in Lasers in Dentistry
(Master Thesis)
2008

RWTH Aachen University
AALZ Aachen Dental Laser Institute

Prof. Dr. Norbert Gutknecht

Table of contents

Acknowledgement

Abstract

1. Introduction

1.1 Diode lasers

1.2 Diode lasers in dentistry

2. Materials and methods

3. Results

3.1 Clinical results

3.1.a 1st case (example for Vision MDL- 10; cw)

3.1.b 2nd case (example for Elexxion Claros; superpulsed)

4. Results of questionnaires

4.1 Patient`s questionnaire

4.1.a Vision MDL- 10 treatment

4.1.b Elexxion Claros treatment

4.2 Surgeon`s questionnaire

4.2.a Vision MDL- 10 treatment

4.2.b Elexxion Claros treatment

5. Summary of results

5.1 Surgeon`s results

5.2 Patient`s results

6. Discussion

7. Conclusion

8. Literature

Acknowledgement

The author wishes to express his special gratitude to Prof. Dr. Norbert Gutknecht, Dr. Rene Franzen and Dr. Jörg Meister for their throughout efforts and patience to guide us on our way to the Master of Science and let us participate in scientific laser dentistry.

Abstract

Purpose :

Dental soft tissue surgery by diode lasers in cw- mode often causes carbonization of the tissues with following necrosis and a delay of wound healing.

In vitro studies have already shown, that superpulsed diode laser surgery has much less disadvantages for the tissues in histological approach.

Purpose of this study is to investigate in vivo, if superpulsed mode of operation can realize an improvement for surgeon and patient in soft tissue surgery.

Materials and methods :

26 patients were treated by diode lasers in different modes of operation for soft tissue surgery. 12 patients were treated by superpulsed Elexxion Claros diode laser :

810 nm; 10- 50 W P_{peak} ; 10- 20 μ s pulse duration; 12000- 20.000 Hz; 400 μ m fiber

and 14 patients were treated by Vision MDL-10 diode laser : 980 nm; 2,5 W; cw- mode and also 400 μ m fiber.

Clinical treatment was documented by photos and questionnaires for patients and surgeons. Questions concerned: carbonization, coagulation, cutting speed, pain, swelling, bleeding, need for drugs, functional reduction and fibrine layer on wounds- during treatment, directly after treatment, after 1 day, after 3 days and after 1 week.

Results :

The clinical observations and the questionnaires showed in most cases significant differences between cw- mode and superpulsed diode laser treatment in surgery.

There was less carbonization in the superpulsed group.

The cutting speed was higher and the cut itself more defined and deeper by using superpulsed mode.

Superpulsed laser treatment had a shorter healing time than cw- mode treatment; the fibrine layer was build faster and also the removal of it was faster.

There was often no swelling after superpulsed diode laser treatment; and if a swelling occurred it was smaller and quicker gone as in cw- mode treatment.

The duration of pain was shorter and the amount of pain smaller in the superpulsed group, therefore the patients in the superpulsed group needed less analgetic drugs.

There was less functional reduction in time and amount for the superpulsed group.

Only coagulation ability was better in the cw- mode group.


Conclusion :

Clinical studies have shown that superpulsed diode laser surgery is superior to continuous wave done treatment. Carbonization and thermal damage of the tissues can be reduced to a minimum, therefore healing is faster as in cw- mode surgery. Generation of a soft tissue cut is faster and more precise. Patients have less pain; in amount and time period. The need of drugs is reduced. There are less functional restrictions and there is less swelling. The advantages of superpulsed mode of operation for soft tissue diode laser surgery are evident. Continuous wave mode should no longer be implemented in diode laser surgery.

1. Introduction

1.1 Diode lasers

Diode lasers are semiconductors, mainly with *GaAs* compounds and mixed crystals. The composition of the mixed crystals, the kind of *GaAs* compound and the temperature determine the wavelength. In diode lasers the laser beam is generated directly at the p-n junction out of electrical energy; no arc- or flashlamps are needed for pumping, therefore the efficiency is about 35- 50 % (76,141).



(Meister J,Franzen R: Diode Lasers, Aachen Institute for Laser Dentistry, M4, ED2006)

There are no real resonators in common sense, the dimensions are particularly small (0,1 dm³/ kw). Because of the slit- diffraction at the p-n junction of the semiconductor the beam profile is elliptical and highly divergent. It has to be collimated before use. The output power from a diode is only in the mW region, so several diodes have to be arranged in arrays, bars and stacks to achieve an output power in the W range. The focussed intensity can reach 10⁴- 10⁶ W/ cm².

Diode lasers for dentistry operate in the near infrared region. The most commonly used wavelengths are 810, 940 and 980 nm (141), because these wavelengths are very well absorbed by pigmented tissues, haemoglobine and melanin, which makes the diodes suitable for soft tissue surgery, endodontics, periodontics and LLLT (57,72,141).

The history of diode lasers started 1962 with the first *GaAs* diode laser, 840 nm, pulsed in fluid *He* or *N₂*. It was developed by IBM, General Electric and Lincoln Laboratories.(155). The first red diode laser was introduced.

1966 the first fiber to guide the laser beam was created.

In 1970 the Bell Laboratories produced the first cw- mode diode laser (155).

In 1995, on the IDS, the first diode laser for dentistry was shown. It was an 810 nm diode laser with cw mode. (133). Peak power was 6 W.

1999 diode lasers with 10- 15 W followed.

2000 a short pulsed diode laser with 810 nm, 20 W, 50 μ s pulse and an average power of 2,0- 6,7 W was introduced, which showed a faster cutting and only 50 % of the necrosis zone of an cw- mode diode laser. (156).

2002 the peak power reached 30 W with a pulse duration of 9 μ s, a frequency of 20000 Hz and a limited average power of 10 W. Today frequencies up to 30000 Hz and peak powers of 50 W are possible. To gain a sufficient result the peak power has to be more than 8 W in this case. (141).


During the first years of diode laser treatment in dentistry only cw mode was possible.

Several studies at that time showed that cw mode and 1 W was enough to reach a bactericidal effect on and in roots and root canals, so as on implant surfaces.

(13- 16,51,88,106,118,120,156). Application of 3- 4 W in cw mode led very fast to carbonization of the soft tissue; the carbonization caused higher absorption followed by a heavy thermal damage and necrosis of the tissue. (71,146,149,150,156).

To approach better results in soft tissue treatment without much carbonization it was necessary to interrupt the cw mode. That was done by chopping the cw mode. Pulses down to several 100 μ s were realized. (141,97,145). The peak power of the pulses was in fact not higher than the peak power of the cw mode pulse, but the applied dose was decreased and the carbonization and thermal damage was reduced.

Operation modes of a diode laser



(Meister J, Franzen R: Diode Lasers, Aachen Institute for Laser Dentistry, M4, ED 2006)

A further progress were pulse durations of $9 \mu\text{s}$, a frequency of 20000 Hz and a peak power of 50 W, which could only be generated by DPL (digital pulse technique).

Using these parameters could generate a fast, sharply edged cut nearly without thermal damage of the surrounding tissue and only minimal carbonization at the cutting edges. (17,24,26,34,35,109,153,154,156).

To make a statement about the laser - tissue interactions something about the general mechanisms of laser irradiation has to be said.

If a laser beam hits the tissue surface 4 different effects can be noticed. (67,68)

1. Reflection
2. Transmission
3. Scattering
4. Absorption

Very important for a surgical effect on the tissue is the absorption of laser light. The kind of tissue and the wavelength are essential for absorption. There are some more factors been involved in laser- tissue interactions, so as the frequency, the dose applied, the pulse duration, power density and very important: the mode of operation (cw, pulsed, superpulsed) which is a main question of the in vivo studies done here.

Dependend on the just mentioned factors there are different effects to be watched in the tissue. (52,67,68,151,152)

Photodynamical effects

Photochemical effects

Thermal effects

(34,64,67,68,71,146,147,148,149,150,151,152)

T= 40- 60 °C: oedema, functional imbalance of metabolism


(up to 50°C reversibel)

T= 60- 100°C: denaturation of the proteins, coagulation, contraction
because of dehydration

T > 100°C : disruption, evaporation of water, ablation

T > 150°C : vaporisation, ablation

T > 300°C : carbonisation



(Franzen R, Meister J: Dosimetry. Aachen Institute for Laser Dentistry, M3, ED 2006)

1.2 Diode lasers in dentistry

As mentioned above diode lasers were used in dentistry since 1995.

Mostly with a wavelength of 810 nm. Later the wavelength 980 nm follows and recently a wavelength of 940 nm is in progress.

In comparison to other lasers the diodes are more economical, the construction is more simple, laser light is generated directly from electrical energy, service periods are rare (maintenance only in intervals of 10000 hours) and the efficiency is high. The laser beam is guided through a quartz- fiber, which makes it easier to reach problematic zones in the oral cavity. Unfortunately diode lasers are only applicable in soft tissue procedures, endodontics, periodontics and LLLT, but not practicable for the ablation of hard tissues. In soft- tissue treatment the diode lasers have similar effects than the Nd:YAG laser, which is not so far away regarding to the wavelength (1064 nm).

Diode lasers are established in endodontical treatment because of the good bactericidal effect (2,13-15,26,28,29,57,72,98).

Diodes are implemented in periodontal treatments and periimplantitis to remove granulated tissue and to sterilize the surface of root and implant (2,11,13,16,18-23,40,50,57,60,72,80,82,85,104,105,134-136,139,140,144).

Also in little surgery diodes are used because of several advantages:

- efficient cutting
- good haemostasis
- nearly without bleeding and therefore good visibility of the op site
- mostly no sutures needed
- no or only small post op oedema
- bactericidal effect
- pain reduction post and intra op treatment
- no secondary bleeding
- precision of cut
- calculated depth of cut
- reduced application of drugs (antibiotics, analgetics..)
- only minimal destruction of adjacent tissues
- uncomplicated handling

- high patient acceptance
- low level of scar forming
- use on patients with haemorrhagic diathesis without or only with little substitution
- good handling due to fiber optics
- reduction of needed instruments (change of parameters- different effects)
- treatment duration shorter (no suture, not often change of instrument..)
- biostimulation of the surrounding tissue

(1,3-10,12,17,25-27,30-39,42-49,57-69,72,74,77,81,86-89,91-97,106,
109-114,116,119-132,138,141,142,145-154,156)

In the beginning soft tissue surgery by diodes was limited because of the low peak powers and the cw- mode as only mode of operation. The doses which were applied to the soft tissue led very fast to carbonization, thermal damage and necrosis of the tissue if the power exceeds 3-4 W in cw- mode. To achieve a better result the dose had to be reduced.

This can be done by decreasing the power, but then there will be no sufficient cutting ability. Another way is to shorten the treatment time, or to split the continuous wave into pulses; and change the frequency of the pulses.

By using mechanical shutters the cw- mode of the diode was changed into a chopped mode, the dose was reduced, but the peak powers of the pulses equaled the peak powers of the cw- mode. Carbonization was reduced for a certain degree but no efficient cutting could be done. The peak power of the pulses was not sufficiently high enough. So the management of the pulses was changed from mechanical to digital. (DPL). Now very short pulses with high peak powers could be generated to avoid carbonization and thermal damage as good as possible. Several studies in the recently past years have shown, that short pulses with high peak powers could generate a much better result than lower powers with longer pulse duration. (17,24,26,34,35,109,153,154)

Already Bach et al. (109) could show 2008 in vitro, that different modes of operation; cw-pulsed- superpulsed; led to different effects.

Result of the study was the suggestion that cw- mode diodes should no longer be used in soft tissue surgery because of the high amount of thermal damage and necrosis of the

adjacent tissue. The superpulsed diode had much less carbonization and thermal damage as side effect. Interestingly 980 nm and 810 nm diodes have been used in that study-- but histologically there was no difference between the 810 nm and 980 nm diode; only the mode of operation was responsible for the grade of cell damage.

Case reports of Maiorano et al. (153) from 2006 could show in vivo that superpulsed diode lasers caused only minimal involvement of the adjacent tissue during surgery, that there was a clean cut without thermal side effects, but still coagulation. Result was an excellent postoperative condition with minimal swelling and pain.

Until now there are unfortunately no clinical or in vivo studies to show the different results of diode laser soft tissue surgery with different modes of operation. Theoretically the difference in treatment results with cw- mode and pulsed is comprehensible; in vitro it is shown already (109) but in vivo it is still unproven. It would be very nice to have no more carbonization during surgical treatment, because carbonization causes a higher grade of absorption and due to that a higher thermal damage of the tissue, which will lead to more discomfort of the patient and a delayed wound healing.

Therefore it is a major task of this thesis to answer the question if highly pulsed diode lasers are more suitable for soft tissue surgery than cw- mode diode lasers and how the achievable results differ.

Of course a number of noticed parameters in the following study, so as patient`s pain, influence on oral hygiene, reduced masticatory or linguallly function etc. are subjective criteria, but at least it can show a tendency for the right direction to go.

2. Materials and methods

Within a period of 5 months 26 dental surgery treatments have been done by diode laser. Predominantly simple laser cuts, but also removal of hyperplasia and fibromas, exposure of teeth and implants, abscess incision and gingivo- and vestibuloplasty.

14 patients have been treated by **Vision MDL-10** diode laser :

- Vision Lasertechnik GmbH, Lügensteinweg 27, 30890 Göxe, Germany



- cw- mode or 20 Hz mode
- max. 2,5 W
- 200 μm and 400 μm fiber
- 980 nm

All procedures have been done with: **cw mode; 2,5 W; 400 μm fiber**

Treatment	P _{peak}	Mode
Single cut (4x)	2,5 W	cw
Incision abscess (1x)	2,5 W	cw
Hyperplasia (1x)	2,5 W	cw
Exposure implant (3x)	2,5 W	cw
Exposure tooth (1x)	2,5 W	cw
Fibroma (2x)	2,5 W	cw
Gingivoplasty (1x)	2,5 W	cw
Vestibuloplasty (1x)	2,5 W	cw

12 patients were treated by ***Elexxion Claros*** diode laser :



- Elexxion AG, Schützenstraße 84, 78315 Radolfzell, Germany
- 810 nm
- 10 mW- 50 W P_{peak}
- 2,5 μs - cw pulse duration
- 200 μm , 300 μm , 400 μm and 600 μm fiber
- digitally superpulsed (8- 20.000 Hz)

Treatment	P _{peak}	Frequency	Pulse duration	P _{avg}
Single cut (3x)	50,0 W	20.000 Hz	11 μ s	11,0 W
Incision abscess (2x)	10,0 W	20.000 Hz	20 μ s	4,0 W
Hyperplasia (1x)	50,0 W	12.000 Hz	10 μ s	6,0 W
Frenectomy (1x)	15,0 W	20.000 Hz	16 μ s	4,8 W
Exposure implants (2x)	15,0 W	15.000 Hz	10 μ s	2,25 W
Exposure tooth (1x)	25,0 W	15.000 Hz	10 μ s	3,75 W
Fibroma (2x)	50,0 W	12.000 Hz	10 μ s	6,0 W

All treatments have been done with the 400 μ m fiber for better comparison.

Both lasers have been used with fibers in contact to tissue.

The fibers had been under permanent control for tidiness, because contaminated fibers decrease the applied energy, the cutting speed drops down, the treatment time prolongs and more thermal energy is given to the tissue (145).

After sufficient anesthesia, for abscess incision only superficial, in all other cases infiltration anesthesia, the treatment was done in consideration of the laser safety directions given for laser treatment of class 4 lasers.

Before, during and directly after laser treatment photos were taken and the patients had to fill in a questionnaire. Another questionnaire had to be done by the surgeon. 1 day, 3 days, 1 week, 2 weeks and sometimes 3 month later there was a recall to control the op sites clinically and to take some more photos. The patients had to fill in their questionnaires until 1 week after surgery. Asked subjects were intra- and postoperative pain, need of analgetic drugs, difficulties in oral hygiene or mastication. The surgeons questionnaire asked about cutting speed of the laser, bleeding intra- and post op, swelling, carbonization, coagulation and fibrine layer covering. Time schedule was the same as for the patients.

Questionnaire

for patient

Name:

Subject : Laser treatment

Pain :

During laser treatment :

0	1	2	3	4	5	6	7	8	9
---	---	---	---	---	---	---	---	---	---

No pain

Unbearable pain

Directly after laser treatment :

0	1	2	3	4	5	6	7	8	9
---	---	---	---	---	---	---	---	---	---

No pain

Unbearable pain

1 day after laser treatment :

0	1	2	3	4	5	6	7	8	9
---	---	---	---	---	---	---	---	---	---

No pain

Unbearable pain

3 days after laser treatment :

0	1	2	3	4	5	6	7	8	9
---	---	---	---	---	---	---	---	---	---

No pain

Unbearable pain

1 week after laser treatment :

0	1	2	3	4	5	6	7	8	9
---	---	---	---	---	---	---	---	---	---

No pain

Unbearable pain

Analgesic drugs needed :

☐

only 1 day

☐

3 days

☐

the first week

☐

more than a week

Reduced masticatory or linguallly function :

☐

only 1 day

☐

3 days

☐

the first week

☐

more than a week

Difficulties in oral hygiene :

☐

only 1 day

☐

3 days

☐

the first week

☐

more than a week

Questionnaire

for surgeon

Name:

Subject : Laser treatment

During laser treatment

Speed of cutting

0

1

2

3

very slow

slow

fast

very fast

Carbonization

0

1

2

3

no carbonization

little carbonization

carbonization

heavy carbonization

Coagulation

0

1

2

3

no coagulation

little coagulation

good coagulation

very good coagulation

Directly post OP

Bleeding

0

1

2

3

No bleeding


little bit

bleeding

heavy bleeding

1 day post OP

Bleeding

			
No bleeding	little bit	bleeding	heavy bleeding

Turgor, swelling

			
No swelling	little swelling	swelling	large swelling

Fibrine layer :

			
No layer	small region	mostly covered	completely covered

3 days post OP

Bleeding :

			
No bleeding	little bit	bleeding	heavy bleeding

Turgor, swelling :



No swelling



little swelling



swelling



large swelling

Fibrine layer :



no layer



small region



mostly covered



completely covered

1 week post OP



bleeding



swelling



fibrine layer

As already mentioned in the introduction the purpose and aim of this in vivo study is to show, that superpulsed diode lasers are more suitable for oral surgery than cw- mode diode lasers, because it is expected that there is nearly no carbonization and thermal damage of the adjacent soft tissue because of the very short pulse durations, the high peak power and high frequency.

Until now there are not enough data from in vivo studies available to support the supposed results, but in vitro studies already came to a similar resume.

(24,34,35,109,132,153)

The diode lasers used in this study are not of the same wavelength; on the one hand there is an 810 nm diode (Elexxion Claros) on the other hand a 980 nm diode (Vision MDL- 10).

According to the wavelength it must be said that normally there is a stronger thermal effect using the 980 nm diode (72) but as shown in a recent in vitro study (109) from 2008, the mode of operation is much more responsible for the clinical results on soft tissue as the diode`s wavelength.

3. Results

3.1 Clinical results

To show the visual results of a soft tissue surgery treatment by diode laser in cw and superpulsed mode 2 similar cases were taken as example for all clinical cases. Treatment procedure was the same in both cases and the local sites were comparable.

In both procedures a fibroma was removed out of the inner lip.

Once on the right and once on the left side of the lower jaw.

For better comparison it would have been perfect if both treatments had been done on the same patient, but there was no such case during the investigation period.

3.1.a 1st case : as example for Vision MDL- 10; cw

Removal of fibroma by Vision MDL 10 diode laser:

980 nm; cw mode; 2,5 W; 400 μ m fiber; contact mode

Female patient, 60 years old, hypertonia. No other general diseases.

18, 15, 25, 28, 36, 48 missing. Bridge 14- 16.

Fair oral care. General horizontal loss of bone. No periodontal problems.

Patient was often sucking on the fibroma.

Fibroma was situated at regio 43, inner lip. Diameter about 8 mm.

After little infiltration anesthesia with Ultracain DS® the fibroma was removed under permanent tension with the 400 μ m fiber in contact.




Before treatment.



During treatment (980 nm; cw; 2,5 W; 400 μ m fiber)

Because of the good absorption of the 980 nm wavelength in melanin and haemoglobine there was a sufficient cutting ability and no bleeding.



Directly after treatment.

The treatment procedure took a little longer because the cutting speed of the 980 nm diode in cw mode with only 2,5 W peak power was slow.

No bleeding because of the very good coagulation of the 980 nm diode, but strong carbonization of the soft tissue was seen.

No sutures needed.

There was no pain or discomfort for the patient during the surgical procedure.

1 day after surgical treatment no secondary bleeding, but a little swelling and a fibrine layer that covered the wound was seen. The patient had no pain but the mastication, the lingual function and the oral hygiene was influenced.

1 day after surgical treatment.



3 days after surgical treatment

No bleeding, a little swelling and a fibrine covered wound was seen.

The patient had no longer functional limitations. No pain.

1 week after treatment
No bleeding and no more fibrine layer was seen.



3 months later. No rezidiv. Total recovery. No scar.

3.1.b 2nd case : as example for Elexxion Claros; superpulsed

Removal of fibroma by Elexxion Claros diode laser.

810 nm; 50 W P_{peak} ; 12000 Hz; 10 μ s pulse duration; 6 W P_{avg} ;
400 μ m fiber; contact mode

Female patient, 72 years old, diabetic disorder. No other general diseases.

18- 28, 38- 35, 32-48 missing.

Total upper denture. Lower partial denture with clasps at 33, 34.

General horizontal and vertical loss of bone. Fair oral care.

Fibroma was situated at regio 35, inner lip, lower jaw. Diameter about 11 mm.

After sufficient anesthesia with Ultracain DS® the fibroma was removed under good tension with the 400 μ m fiber in contact.



Before treatment (little bleeding spot from injection of anesthesia)

During treatment

(810 nm; 50 W P_{peak} ; 12000 Hz; 10 μs pulse; 6 W P_{avg} ; 400 μm fiber)



No bleeding because of good coagulation of 810 nm wavelength (good absorption in haemoglobine and melanin).

Extremly fast cutting speed because of high power, short pulses and high frequency.

The patient felt no discomfort or pain during the surgical treatment time.

Directly after treatment no bleeding of the wound, no sutures needed.

Good coagulation and only a little carbonization zone was seen.

Directly after treatment



1 day after treatment

3 days after treatment

No bleeding, no swelling. Wound is completely covered by fibrine layer.

Patient had no discomfort or pain.




1 week after treatment

No bleeding, no swelling, partial fibrine layer. Patient without any discomfort.

2 weeks after treatment

Nearly full recovery. Only a little impression of the soft tissue remained.

No rezidiv, no scar. Healing much quicker than with cw diode laser.



4. Results of questionnaires

In total there were 52 questionnaires to analyse.

26 from the patients and 26 from the surgeons.

4.1 Patient`s questionnaire :

4.1.a VISION MDL- 10 treatment ➡ 14 patients; cw

During laser treatment

All patients felt no **pain** during laser treatment (100 %)

Directly after laser treatment

All patients without any **pain** (100 %)

1 day after laser treatment

Number of patients	Pain: scale 0- 9	Percentage of total number of patients
3	0 (no pain)	21,43%
3	1	21,43%
5	2	35,71%
2	3	14,29%
1	7 (strong pain)	7,14%

3 days after laser treatment

Number of patients	Pain: scale 0- 9	Percentage of total number of patients
10	0 (no pain)	71,43%
3	1	21,43%
1	7 (strong pain)	7,14%

1 week after laser treatment

Number of patients	Pain: scale 0- 9	Percentage of total number of patients
13	0	92,86%
1	5	7,14%

Analgesic drugs needed

Number of patients	Analgesic drugs needed	Percentage of total number of patients
8	No	57,14%
5	Yes, 1 day	35,71%
0	Yes, 3 days	0%
1	Yes, 1 week	7,14%

Reduced masticatory or linguallly function

Number of patients	Time of reduction	Percentage of total number of patients
0	No reduction	0%
4	1 day	28,57%
8	3 days	57,14%
1	1 week	7,14%
1	more than a week	7,14%

Difficulties in oral hygiene

Number of patients	Time of handicap	Percentage of total number of patients
0	No	0%
5	1 day	35,71%
8	3 days	57,14%
0	1 week	0%
1	more than a week	7,14%

4.1.b ELEXION CLAROS treatment ➔ 12 patients; superpulsed

During laser treatment

All patients felt no **pain** during laser treatment. (100 %)

Directly after laser treatment

All patients without any **pain** (100 %)

1 day after laser treatment

Number of patients	Pain: scale 0- 9	Percentage of total number of patients
5	0 (no pain)	41,67%
3	1	25%
3	2	25%
1	3	8,33%

3 days after laser treatment

Number of patients	Pain: scale 0- 9	Percentage of total number of patients
8	0 (no pain)	66,67%
4	1	33,33%

1 week after laser treatment

Number of patients	Pain: scale 0- 9	Percentage of total number of patients
12	0 (no pain)	100%

Analgesic drugs needed

Number of patients	Analgesic drugs needed	Percentage of total number of patients
8	No	66,67%
3	Yes, 1 day	25%
1	Yes, 3 days	8,33%

Reduced masticatory or linguall function

Number of patients	Time of reduction	Percentage of total number of patients
2	No reduction	16,67%
4	1 day	33,33%
6	3 days	50%

Difficulties in oral hygiene

Number of patients	Time of handicap	Percentage of total number of patients
1	No	8,33%
5	1 day	41,67%
6	3 days	50%

4.2. Surgeon`s questionnaire

4.2.a VISION MDL- 10 treatment ➡ 14 patients; cw

During laser treatment

Speed of cutting was all time low (100 %)

Carbonization

Number of patients	Degree of carbonization	Percentage of total number of patients
0	no carbonization	0%
2	little carbonization	14,29%
7	carbonization	50%
5	heavy carbonization	35,71%

Coagulation

Number of patients	Degree of coagulation	Percentage of total number of patients
0	no coagulation	0%
0	little coagulation	0%
8	good coagulation	57,14%
6	very good coagulation	42,86%

Directly post op

Number of patients	Degree of bleeding	Percentage of total number of patients
11	no bleeding	78,57%
3	little bleeding	21,43%

1 day post op

There was no recognizable **bleeding** (100 %)

Swelling

Number of patients	Degree of swelling	Percentage of total number of patients
4	no swelling	28,57%
9	little swelling	64,29%
1	swelling	7,14%

Fibrine layer

Number of patients	Size of fibrine layer	Percentage of total number of patients
4	no layer	28,57%
5	small region	35,71%
5	mostly covered	35,71%

3 days post op

There was no recognizable **bleeding** (100 %)

Swelling

Number of patients	Degree of swelling	Percentage of total number of patients
8	no swelling	57,14%
6	little swelling	42,86%

Fibrine layer

Number of patients	Size of fibrine layer	Percentage of total number of patients
2	no layer	14,29%
3	small region	21,43%
3	mostly covered	21,43%
6	completely covered	42,86%

1 week post op

No **bleeding** (100 %) left; no **swelling** (100 %) left.

6 patients of 14 still had a little **fibrine layer** (42,86 %).

4.2.b ELEXION CLAROS treatment ➔ 12 patients; superpulsed

During laser treatment

Speed of cutting

Number of patients	Speed of cutting	Percentage of total number of patients
0	slow	0%
3	fast	25%
9	very fast	75%

Carbonization

Number of patients	Degree of carbonization	Percentage of total number of patients
4	no carbonization	33,33%
4	little carbonization	33,33%
4	carbonization	33,33%
0	heavy carbonization	0%

Coagulation

Number of patients	Degree of coagulation	Percentage of total number of patients
0	no coagulation	0%
2	little coagulation	16,67%
10	good coagulation	83,33%
0	very good coagulation	0%

Directly post op

Bleeding

Number of patients	Degree of bleeding	Percentage of total number of patients
7	no bleeding	58,33%
5	little bleeding	41,67%

1 day post op

There was no recognizable **bleeding** (100 %)

Swelling

Number of patients	Degree of swelling	Percentage of total number of patients
5	no swelling	41,67%
7	little swelling	58,33%
0	swelling	0%

Fibrine layer

Number of patients	Size of fibrine layer	Percentage of total number of patients
2	no layer	16,67%
6	small region	50%
4	mostly covered	33,33%

3 days post op

There was no recognizable **bleeding**. (100 %)

Swelling

Number of patients	Degree of swelling	Percentage of total number of patients
11	no swelling	91,67%
1	little swelling	8,33%

Fibrine layer

Number of patients	Size of fibrine layer	Percentage of total number of patients
0	no layer	0%
4	small region	33,33%
5	mostly covered	41,67%
3	completely covered	25%

1 week post op

No **bleeding** (100 %); no **swelling** (100 %) left.

4 of 12 patients still had a little **fibrine layer** (33,33 %)

5. Summary of results

5.1 Surgeon`s results

The visible clinical results could be described as followed :

During surgical treatment and directly after treatment the soft tissues at the involved area were more carbonized by the Vision MDL- 10 laser, which ran in cw- mode;
35,71 % heavy carbonization.

The superpulsed laser , Elexxion Claros, showed sometimes carbonization too, but the amount of carbonization was much smaller, the color of the carbonized zones was more brown than black and in a third of all cases there was no carbonization at all.

Carbonization	cw- mode	superpulsed
no	0%	33,33%
little	14,29%	33,33%
medium	50%	33,33%
heavy	35,71%	0%

Coagulation was good with both lasers, in some situations (a patient with anticoagulant therapy - Marcumar® - INR > 2,5 ; not substituted) the 980 nm diode in cw- mode (Vision MDL- 10) had a better coagulation.

Coagulation	cw- mode	superpulsed
little	0%	16,67%
good	57,14%	83,33%
very good	42,86%	0%

Both lasers did not allow common bleedings post op, but in some cases there was a little oozing bleeding; more after superpulsed laser therapy than after cw- mode treatment.

Bleeding	cw- mode	superpulsed
post op : no	78,57%	58,33%
post op : little	21,43%	41,67%


The **speed of cutting** was very different, the superpulsed Elexxion diode laser cutted much easier and faster than the Vision diode laser.

The surgical procedure was done much faster by using the superpulsed laser.

Speed of cutting	cw- mode	superpulsed
low	100%	0%
fast	0%	25%
very fast	0%	75%

On pictures taken after the surgical procedure the wound areas looked much nicer when the treatment was done by superpulsed Elexxion laser.

The cutting margins had sharper edges, the cut was more straight and seemed to go deeper into the tissue with one movement.



cw mode



superpulsed

The following days there was a swelling of the treated soft tissue in different sizes noticed for both laser systems.

Swelling 1 day post op	cw- mode	superpulsed
no swelling	28,57%	41,67%
little swelling	64,29%	58,33%
swelling	7,14%	0%

The swelling tendency of the superpulsed diode laser was lower and there was no swelling in significant more cases.

Swelling 3 days post op	cw- mode	superpulsed
no swelling	57,14%	91,67%
little swelling	42,86%	8,33%

After 3 days there was only in 8,33 % a little swelling left after superpulsed diode laser treatment but in the cw- mode group 42,86 % of the patients still had a recognizable swelling. One week after laser treatment there was no swelling left in both laser groups.

Another phenomenon to be watched was the fibrine layer.

Fibrine layer 1 day post op	cw- mode	superpulsed
no layer	28,57%	16,67%
small region	35,71%	50%
mostly covered	35,71%	33,33%

Including 4 osteotomy cuts for cw- mode and 3 for superpulsed mode.

In total there was more fibrine layer after superpulsed treatment.

Fibrine layer 3 days post op	cw- mode	superpulsed
no layer	14,29%	0%
small region	21,43%	33,33%
mostly covered	21,43%	41,67%
completely covered	42,86%	25%

3 days post op the fibrine covering of the wound areas had increased in both groups.

Fibrine layer 1 week post op	cw- mode	superpulsed
no more layer	57,14%	66,67%
still little layer	42,86%	33,33%

1 week post op there was more reduction of fibrine layer in the superpulsed group, while the cw- mode group still had a partial fibrine covered area in 42,68 % of the patients.

5.2 Patient`s results

During laser treatment and directly after laser treatment no patient felt pain or discomfort because there was anesthesia given in every treatment case.

Pain 1 day post op: 0= no pain 9= unbearable	cw- mode	superpulsed
0	21,43%	41,67%
1	21,43%	25%
2	35,71%	25%
3	14,29%	8,33%
7	7,14%	0%

There were more patients without pain and the pain strength was on a lower level in the superpulsed group.

Pain 3 days post op: 0= no pain 9= unbearable pain	cw- mode	superpulsed
0	71,43%	66,67%
1	21,43%	33,33%
7	7,14%	0%

Now there were more patients without any pain in the cw- mode group and also patients with very mild pain were less in the cw- mode group. But still there was a special patient with continued heavy pain.

Pain 1 week post op: 0= no pain 9= unbearable pain	cw- mode	superpulsed
0	92,86%	100%
5	7,14%	0%

After 1 week all patients in both groups were free of pain, the only exception was the special patient of the cw- mode group.

For better understanding it has to be said that the mentioned patient was very algesic; a vestibuloplasty had been done and the extended denture was incorporated directly after laser treatment, so the healing was inhibited for a certain degree and the denture was grinding on the wound area all the time.

Analgesic drugs needed	cw- mode	superpulsed
No	57,14%	66,67%
Yes, 1 day	35,71%	25%
Yes, 3 days	0%	8,33%
Yes, 1 week	7,14%	0%

The need for analgesic drugs was higher in the cw- mode group, 42,85 %.

Need for medication in the superpulsed group: 33,33 % .

The above mentioned special patient needed a week of medication.

Reduced masticatory or linguallly function.

Reduced function	cw- mode	superpulsed
No	0%	16,67%
1 day	28,57%	33,33%
3 days	57,14%	50%
1 week	7,14%	0%
more than 1 week	7,14%	0%

The influence of the laser treatment on functional aspects as e.g. mastication lasted in average 3 days.

In the cw- mode group even longer in 2 cases.

Only in the superpulsed group 16,67 % were without any functional limitations.

Difficulties in oral hygiene	cw- mode	superpulsed
No	0%	8,33%
1 day	35,71%	41,67%
3 days	57,14%	50%
1 week	0%	0%
more than 1 week	7,14%	0%

Nearly the same result for difficulties in oral hygiene.

Average was as well 3 days.

Even here there were patients without difficulties in the superpulsed group; 8,33 %.

6. Discussion

Classical surgical soft tissue treatment by scalpel has its relevancy, but in the last years soft tissue surgery by diode laser is in progress.

There are a lot of advantages for the surgeon:

- efficient and precise cutting with a calculated depth of cut
- good haemostasis- so nearly or completely without bleeding- and therefore better visibility of the place of interest.
- mostly no sutures needed
- only minimal destruction of the adjacent tissue (in pulsed mode)
- uncomplicated handling due to fiber and variety of parameters
- shorter treatment time.

The patients have a lot of advantages also:

- they mustn't be sutured
- no or only little post op oedema
- bactericidal reduction of the wound area
- no secondary bleeding
- low level of scar forming
- patients with haemorrhagic diathesis can be treated without or with only little substitution
- shorter treatment time
- biostimulation of the surrounding area
- reduced application of drugs because of pain reduction intra- and post op
(1,3-10,12,17,25-27,42,43,57,72,106,129,141,145,146,153).

The reduction of pain during laser treatment and the better acceptance of laser treatment in comparison to conventional treatment could be proofed yet for restorative dentistry (Matsumoto 1996, Keller- Hibst 1997-1998, Evans 2000, Takamori 2003) but in this cases the used lasers were not diode lasers.

Discussing the results from the study one certain thing has to be remembered.

The cw mode laser (Vision MDL- 10) had always been used with 2,5 W and cw- mode.

So peak- and average power had always be the same; in every procedure.

That`s the main difference to the superpulsed laser (Elexxion Claros) which had been used with different parameters ; 10- 50 W peak power; 2,25- 11,00 W average power ; 12000- 20000 Hz; 10- 20 μ s pulse duration.

The precise treatment parameters of the Elexxion Claros laser were :

P_{peak}	P_{avg}	frequency	pulse duration
6x 50 W	3x 11,00 W	6x 20000 Hz	6x 10 μ s
1x 25 W	3x 6,00 W	3x 15000 Hz	3x 11 μ s
3x 15 W	1x 4,80 W	3x 12000 Hz	1x 16 μ s
2x 10 W	2x 4,00 W		2x 20 μ s
	1x 3,75 W		
	2x 2,25 W		
Ø 32,5 W	Ø 6,00 W	Ø 16750 Hz	Ø 12,42 μ s

The used peak power was sometimes 20 times higher than the peak power of the cw mode laser; and even the average power was more than double (2,4 times higher). Regarding to this there was no doubt that the Elexxion laser would cut much easier and faster than the Vision laser. On the other hand the used powers were so high that normally a cw mode diode laser user would expect terrible thermal damage and necrosis of the tissue with a lot of carbonization (71,146,149,150,156). But thanks to the superpulsed mode of operation there was less carbonization and less damage of the adjacent tissue, which had been proven histologically in vitro by Bach et al. in 2008 (109). Bach used a cw mode diode laser, a pulsed mode diode laser and a superpulsed diode laser with P_{peak} 30 W; P_{avg} 9,99 W; 20000 Hz to make soft tissue cuts in pig jaws . Unfortunately he didn`t mention the pulse duration. Peak power was nearly the same as in this study, average power 1,67 times higher as in this study and the average frequency also higher. Still the results were much better for the superpulsed laser than for the cw- mode laser.

Maiorana and Salina (153) did a clinical study with a superpulsed laser on 3 patients in 2006. Parameters were: 23,33 W peak power in average; 16667 Hz in average; 10 μ s pulse duration in average; P_{avg} not mentioned.

The conclusion was : clean cut without thermal side effect; instant coagulation; excellent post op conditions with minimal pain and swelling; minimal involvement of the adjacent tissues during surgery.

The parameters were similar to this study, only the peak power was significantly lower and the number of patients was only 3. But the summary of both studies correspond with the results achieved in this study.


The results out of this comparative in vivo study are sometimes not very different for the used lasers, but sometimes very significant, so in case of the

cutting speed :

The use of high powers, short pulse durations and high frequencies offered a high speed for cutting and a deeper cut as it did in cw-mode. The margins of the cut were more defined and more straight using superpulsed mode. This was already described several times in literature. (26,34,145). When using the superpulsed diode with high power and frequency for cutting the first time, the normal treatment behaviour had to change; the movement speed of the fiber had to be increased because of the rapid cutting velocity. When using the cw- diode laser the cut was flattish and the cutting speed was low. Certainly the low peak power of the Vision MDL- 10 was one great inhibition factor for speed. Another one may have been the wavelength of 980 nm. As known the cutting efficiency of an 810 nm diode laser is better than it is with a 980 nm diode

(Gutknecht N; script M6- Master of Science in Lasers in Dentistry,ED2006; AALZ- Aachen; 02/2008).

The 810 nm diode has less penetration into the depth, a lower absorption in water, a lower absorption in HbO₂ and nearly the same absorption in Hb as a 980 nm diode. The 980 nm diode creates more thermal energy at the surface and penetrates more into depth. The zone of necrosis is larger and it is more dangerous for cutting.



(Pict. from Ezlase brochure; Biolase Technologie Inc. 2007)

Surprisingly Bach et al. (109) could prove 2008 in in vitro studies, that there was no difference from the histological point of view between 810 nm and 980 nm diode lasers, when both are used in the same mode of operation.

If we assume that Bach`s results were correct and we neglect the difference between in vitro and in vivo studies, we can conclude together with the attained results from this study that superpulsed mode of operation will lead to faster cutting speed and a more precise cut with less tissue damage.

There are a lot more parameters apart from the wavelength, power, frequency, pulse duration, fiber diameter and mode of operation that will influence the cutting ability,

so as e.g. : kind of tissue, pigmentation, race, blood circulation in tissue, applied dose, treatment time etc. but these parameters are of subsidiary relevance in this study.

The above shown reference cases are done with patients of the same race and same gender. The soft tissue region was the same, both patients with similar pigmentation. The treated sites were mirror- inverted at the inner lip of the lower jaw.

Carbonization :

Next point to view is the carbonization of the tissue. If there is a lot of carbonization the destruction of the surrounding tissue is large (37,45,71,109).

Carbonization is changing the absorption of the treated tissue; it is increasing due to the dark color which absorbs the diode laser light much better than light colors do.

This means there is more thermal energy applied in the surface and necrosis of the adjacent structures is created.

As we can see from the study the changing of parameters could influence the degree of carbonization and coagulation; which had already been described in literature(145).

In 1999 Goharkhay et al.(7) stated that the horizontal and vertical dimension of the tissue destruction is neither addicted to the diameter of the fiber, nor related to the mode of operation (cw or pulsed) but is only depending on the average power used.

Both assumptions could be rebut until today. Using the same power settings the dose applied in a defined treatment time is lower with a larger fiber diameter. In using smaller fiber diameters the applied energy is distributed in a smaller area, the dose is higher, the effect is larger, there is more thermal energy brought into the tissue.

The risk of carbonization is increasing.

So the tissue destruction is certainly related to the diameter of the fiber.

The used average power in case of the Vision MDL-10 laser was 2,5 W; in case of the Elexxion Claros 6,0 W.

According to Goharkhay the tissue damage must be larger when using the Elexxion Claros

because of the higher average power, but it was just the opposite, because the mode of operation was changed.

Also in literature a much smaller degree of carbonization and destruction is described when using pulsed instead of cw- mode (24,26,34,35,41,61,62,109,145,153).

This can be consolidated by the results achieved here.

The cw- mode laser produced carbonization in every treatment case, 100 % .

The superpulsed diode laser only in 66,66 % of the cases.

No patient with heavy carbonization in the superpulsed group.

In the cw- mode group, 35,71 % had severe carbonization,

The total amount of carbonization was much higher in the cw- mode group.

Maybe carbonization could have been less in the superpulsed group regarding to the treatment practice, because the fiber had to be moved much faster as in cw- mode and the surgeon had to adapt to the new procedure technique.

The fiber diameter is another factor for the applied dose and perhaps the carbonization could have been decreased by using a larger fiber; but for better comparison both lasers were used with 400 μm fibers. All other parameters which could influence the degree of carbonization, such as peak power, average power, frequency, pulse duration etc. were accepted in the way the manufacturers had programmed the diode laser units.

Exception was, as just said, only the fiber diameter of 400 μm which was prescribed by the manufacturer of Vision MDL-10 diode laser for surgery. In case of the Elexxion Claros the prescribed fiber diameters should have been changed from 200 to 400 or 600 μm according to the planned surgical treatment.

Coagulation :

Next aspect to discuss is the coagulation ability of both modes of operation.

In literature the coagulation efficiency is described as good (17,30,42,67,68,131,145) to very good (7,37,51,72,142) for cw- mode.

For superpulsed mode the coagulation mentioned is very good (17,35,143) and good with minimal bleeding (34).

Geldi C et al. (35) described 2006 in a study, that superpulsed diode lasers have a larger coagulation area but smaller carbonization area as cw- mode diode lasers.

In this study the gained results from the surgeons questionnaire demonstrated in direct comparison that coagulation was stronger by using cw- mode.

During treatment the coagulation was very good in 42,86 % of the treated patients and good in 57,14 %. For superpulsed there was no very good coagulation, but 83,33 % good and 16,67 % little coagulation which meant that a few patients had a minimal bleeding. In cw- mode there was none.

In literature often no bleeding after surgical treatment is mentioned (17,26,29-31,37,72) but in this study there was an oozing bleeding in 21,43 % of the cw- group patients after surgery and also in 41,67 % of the superpulsed group after finishing the surgical treatment , which is a significant difference.

This might have been a derivation of the higher thermal damage of the cw- mode laser. The blood vessels were sealed better by the influence of the larger amount of thermal energy delivered by the cw- mode of operation.

Therefore a patient with anticoagulant therapy (Quick 20 %; INR > 2,5) was treated with the cw- mode diode laser without any substitution; and there was no bleeding at all.

Mentionable in comparison to other studies is that 7 of 26 patients have been treated with a simple diode laser cut for a following osteotomy.

These patients were sutured after osteotomy and therefore a minimal oozing bleeding occurred. This has to be kept in mind for better understanding of the given results.

Swelling :

Swelling of the wound area and adjacent tissue 1 day after surgery was significantly less in the superpulsed group, 41,67 % with no swelling.

Only 28,57 % with no swelling in the cw- group.

The size of swelling didn't differ much in both groups.

3 days after treatment there were :

Already 91,67 % of the patients free of swelling in the superpulsed group.

Only 57,14 % of the cw-mode group; which is a significant difference.

The reason for this could be the degree of tissue damage caused by laser irradiation.

Regarding to the lower degree of carbonization and thermal damage in the superpulsed group it was explainable that there was less swelling.

The literature review showed different kinds of results which varied from no swelling (17,72) to little swelling (26,37,41,46,142,145,153) for both kinds of lasers.

Swelling is a result of the inflammation reaction of the body and, often described in literature, the inflammation reaction is more serious using diode laser surgery than it is by using a traditional scalpel (26,31,32).

Pain :

A second characteristic sign for inflammation is pain.

During and directly after laser treatment there was no pain because all patients had been given local anesthesia before surgery. Surgical treatment with a diode laser always requires anesthesia, because the thermal energy applied always generates pain.

In literature classifications as mild pain (37), less pain (26,34,145), reduced pain (30,41,46,72,142,145), minimal pain (153) and no pain (17) are found and all authors agree that the degree of pain is reduced after laser treatment in comparison to classical treatment by scalpel.

Looking to the results of the questionnaire a recognizable larger percentage of patients had no pain 1 day post op in the superpulsed group: 41,67 % .

Only 21,43 % in the cw group had no pain.

Nearly double of the patients in the superpulsed group were painfree.

The remaining 58,33 % of the superpulsed group only had pain on the 1st - 3rd niveau.

In the cw group 78,57 % of the patients were also on the 1st - 3rd niveau.

Only 1 patient (the very algesic one, as described above) was on the 7th niveau of pain.

After 3 days the situation was completely different.

Now more patients without pain could be found in the cw group, 71,43 %, while only 66,67 % of the superpulsed group were painfree.

But the difference of the percentages was not as significant as it was after the first day.

The rest of the patients could be found on pain niveau 1 in both groups, exception, of course, the very algesic patient, still on pain niveau 7.

Saaleh et al. (37) noticed an average of 3 days of pain for all patients in their study after cw- mode treatment.

1 week after surgical laser treatment no patient felt pain any longer (exception: the very algesic patient on niveau 5).

Need for drugs :

After laser treatment the patients had been supplied by an analgesic drug, Ibuprofen 400 mg.

42,85 % of the patients in the cw group needed medication, but only 33,33 % of the patients in the superpulsed group.

With a high percentage only 1 day.

That corresponded to the findings for pain that have been done before.

In general we can say that the pain sensation was less in the superpulsed group.

These results and the following conclusions for reduced masticatory and linguallly function so as the difficulties in oral hygiene were subjective and will differ by treating other patients because the degree of sensation is very different from person to person.

There will be no reproducable results but only a hint to the right direction.

Reduced masticatory and linguallly function :

In order to the traumata set during laser treatment there was a functional reduction for masticatory and/ or linguallly function in both groups.

In the cw group every patient was affected.

In the superpulsed group at least 16,67 % of the patients had no functional limitations.

50- 57,14 % of the patients were handicapped for 3 days in both groups.

That corresponded once again to the study of Saaleh et al. (37), who described a 3 day period of painful mastication and speech after cw- mode treatment.

Difficulties in oral hygiene :

Nearly the same results for oral hygiene.

In average 3 days of limitation, but this time

only 8,33 % of the patients in the superpulsed group without any difficulties.

In the cw group all patients with difficulties.

Now three of the five signs for inflammation were reviewed. (dolor, tumor, functio laesa)

2 others, calor and rubor were not registered during the study.

But a sign of healing, better secondary healing, should be discussed subsequently.

The Fibrine layer :

The building of fibrine is a part of the secondary haemostasis, a part of the plasmatic haemostasis. The thrombocytes were connected by a dense network of fibrinous fibers.

The grade of fibrine layer relies in a certain extent to the grade of healing.

The faster the fibrine layer is removed, the faster the wound healing is in progress.

2- 4 hours after trauma the extravasation starts , a coagulum is build

with blood- and plasmacells on the wound surface. Fibrin fibers are connecting the thrombocytes and build a dense network of fibers. This normally lasts to the 4th day .

Then proliferation starts, granulated tissue is build and the fibrin fibers are reduced.

From the 5th day on regeneration starts and there is no more fibrine layer.

The faster the fibrine layer is removed, the faster the wound healing is in progress.

The observed results regarding the fibrin layer differed from the normally seen wound healing after classical treatment by scalpel.

In this study 26 patients were treated, but 7 of them only by simple laser cut for osteotomy. The cut was sutured afterwards, so there was primary wound healing and no visible fibrine layer could be seen.

The first day no fibrine layer was seen on the wounds in 28,57 % of the cw group and only 16,67 % of the superpulsed group.

A small region was covered in 50 % of the superpulsed, but only in 35,71 % of the cw group.

Mostly covered wounds were noticed in 33,33 % (superpulsed) and 35,71 % (cw) in both groups.

This meant that in the superpulsed group already 83,33 % of the wounds had a fractional fibrine layer. In the cw group only 71,42 % of the wounds were partially covered with fibrine.

After 3 days all wounds; 100 % were covered with fibrine in the superpulsed group: 41,67 % were mostly covered; 33,33 % partially and 25 % completely.

In the cw group already 42,86 % were completely covered; 21,43 % mostly; 21,43 % partially, but still 14,29 % were not covered yet.

So far fibrine layer building was faster and more complete in the superpulsed group.

In comparison to healing after conventional surgical treatment there was a delay of healing time, because fibrine layer building started later and needed more time to cover the wounds. In consequence of delayed wound healing there was still a fibrine layer seen after 1 week in 33,33 % of the superpulsed group and in 42,86 % of the cw group.

Delayed wound healing after surgical diode laser treatment has been described in literature already (26,31,32) and as we could see from the results above the wound healing of the cw group started later and lasted even longer as the prolonged wound

healing of the superpulsed group.

This might have been an effect of the stronger thermal damage and increased carbonization in the cw group. The cells and blood vessels in the adjacent tissue were much more damaged and needed longer to reorganize for healing.

Superpulsed treatment left more intact tissue at the wound margins, so healing could begin earlier as in cw mode.

Even if there was an improvement in healing time there was still a difference to normal healing, where fibrine layer removal begins already at the end of the 4th day.

So the intactness of the tissue seemed to be essential for regular wound healing.

Superpulsed diode laser treatment did not cause much carbonization, but there was as well coagulation, which meant applied thermal energy on the blood vessels and cells; causing damage and therefore delayed healing.

As option for surgical laser treatment an Er:YAG laser with waterspray and air can be used instead. The Er:YAG laser also generates a soft tissue cut, but with a resulting bleeding. There will be no carbonization and nearly no coagulation, therefore no or only minimal damage of the tissue and as effect a much quicker and nicer wound healing and less discomfort for the patient as after diode laser treatment.

In summary we can say that superpulsed diode laser treatment is contrary to cw- mode diode laser treatment much more gentle for the tissue (less damage, quicker healing), more comfortable for the patient (less pain and swelling) and an improvement for the surgeon (more precise and deeper cut, less burnt smell, shorter treatment time).

Only the coagulation ability of the superpulsed diode laser seems to be minimal inferior of the cw- mode diode laser.

In combination with the findings of the Bach study (109) and the case reports of Maiorana and Salina (153) we can conclude that superpulsed diode lasers should be the state- of- art treatment nowadays if there is spoken about surgical diode laser treatment.

Of course these are only 2 further studies and the results have to be reconfirmed by various studies with a larger amount of patients, but the direction to go is clearly given.

7. Conclusion

Clinical in vitro and in vivo studies have shown that soft tissue surgery done by superpulsed diode lasers is more reasonable as doing the same procedures by cw- mode diode lasers. Carbonization and thermal damage of the adjacent tissue can be reduced to a minimum, the soft tissue cut can be generated faster, the cut is more precise.

Because of the lower amount of tissue destruction the healing is faster as in cw- mode.

It is more comfortable for the patients regarding to the post op pain and swelling, the patients will need less drugs and the functional abilities are not reduced as much as in continuous wave mode.

In terms of Bach`s histological investigation of the treated soft tissue (109) we probably have to share his opinion, that cw- mode diode lasers are no longer suitable for soft tissue cutting in dentistry.

The advantages of the superpulsed diode lasers in soft tissue surgery are evident and with the meanwhile achieved peak powers, the very short pulse durations and the high frequency the surgical treatment can be improved and probably the range of treatment can be expanded.

8. Literature

001. Gouw- Soares S, Tanji E, Haypek P, Cardoso W, Eduardo CP:
The use of Er:YAG, Nd:YAG, and Ga-Al-As-Lasers in Periapical Surgery:
a three year clinical study. J Clin Lasers Med Surg 19(4): 193-198,2001
002. Schoop U, Kluger W, Moritz A, Sperr W: The bactericidal Effect of Different
Laser Systems in Deep Dentin Layers. Laser Congress, Florence 2003
003. Wigdor HA, Walsh JT, Featherstone JDB, Visuri SR, Fried D, Waldvogel JL:
Lasers in Dentistry. Laser Surg Med 16: 103-133, 1995
004. Pick RM, Colvard MD: Current status of Lasers in soft tissue dental surgery.
J Periodontol 64(7): 589-602, 1993
005. Horch HH: Laser in der oralen Chirurgie.
Vortrag beim 11. Internationalen Jahreskongress der DGL. Berlin 2002
006. Luomanen M, Virtanen I: Healing of laser and scalpel incision wounds of rat
tongue mucosa as studied with cytokeratin antibodies.
J Oral Pathol 16: 139-144, 1987
007. Goharkhay K et al.: Effects on oral soft tissue produced by a Diode laser in
vitro. Laser Surg Med 25: 401-406, 1999
008. Goharkhay K, Moritz A, Schoop U, Pattera C, Rumetzhofer A, Wernisch J,
Sperr W: Auswirkung unterschiedlicher Laserwellenlängen auf die orale
Schleimhaut- eine inVitro Studie. Stomatol 97: 173-179, 2000
009. Stock K, Hibst R: Comparative in vitro study of Diode Laser (940 nm) and
Nd:YAG laser for oral surgery. Lecture Laser Congress, Florence 2003
010. Rosenberg C, Tadir Y et al.: Endometrial laser ablation in rabbits : a comparative
study on three laser types. Laser Surg Med 10: 66-73, 1990
011. Romanos G et al.: Laseranwendung in der Parodontologie- aktueller Stand.
Parodontologie 4: 299-312, 1998
012. Gaspar L : The use of high power lasers in oral surgery.
J Clin Laser Med Surg 12(5): 281-285, 1994

013. Moritz A, Gutknecht N, Dörtbudak O, Schoop U, Schauer P, Sperr W: Bacterial reduction in periodontal pockets through irradiation with a Diode laser: a pilot study. *J Clin Laser Med Surg* 15(1): 33-37, 1997
014. Moritz A, Gutknecht N, Goharkhay K, Schoop U, Sperr W: In-vitro Irradiation of infected Root Canals with a Diode Laser: Results of Microbiologic, Infrared Spectrometric and Stain Penetration Examination. *Quintessence Int* 28(3): 205-209, 1997
015. Gutknecht N, Moritz A, Conrads G, Lampert F: The Diode Laser and its bactericidal effect in root canals: an in vitro study. *Endodontie* 3: 217-222, 1997
016. Hartmann HJ, Bach G: Diodenlaser- Oberflächen- Dekontamination in der Periimplantitistherapie. Eine Drei- Jahres- Studie. *ZWR* 106: 524-526, 1997
017. Romanos G, Nentwig GH: Diode Laser (980 nm) in oral and maxillofacial surgical procedures: observations based on clinical applications. *J Clin Laser Med Surg* 17(5): 193-197, 1999
018. Romanos GE: Treatment of periimplant lesions using different laser systems. *J Oral Laser Appl* 2: 75-81, 2002
019. Romanos GE, Everts H, Nentwig GH: Effects of the diode (980 nm) and Nd:YAG (1064 nm) laser irradiation on titanium discs. A SEM examination. *J Periodontol* 71: 810-815, 2000
020. Salina S, Beretta M, Speroni S, Maiorana C: Treatment of peri- implantitis with Diode Laser and GBR: preliminary results.lecturs. *Laser Congress, Florence* 2003
021. McDavid VG, Cobb CM, Rapley JW, Glaros AG, Spencer P: Change in temperature of subjacent bone during soft tissue laser ablation. *J Periodontol* 69: 1278-1282, 1998
022. Soibelman M: Verwendung dentaler Lasersysteme in der zahnärztlichen Praxis. *Laser Journal* (4): 30-33, 2007
023. Moritz A, Gutknecht N, et al.: Keimreduzierung in Zahnfleischtaschen durch Bestrahlung mit einem Diodenlaser. *J Clin Laser Med Surg* 15(1): 33-37, 1997

024. Neckel C: Eine klinische und histologische Untersuchung zum Thema CW-Mode oder 10.000- Hertz- topgepulster Ga-Al-As- Diodenlaser. DZW 9, 2001
025. Bach G, Gutknecht N, Schneider H: Lasertherapie in der Zahnheilkunde. Diodenlaser in der Praxis. Apollonia Verlag, Linnich, 1998
026. Moritz A, et al.: Orale Lasertherapie. Quintessence Verlag, Berlin 2006
027. Frentzen M, Koort HJ: Lasertechnik in der Zahnheilkunde Dtsch Zahnärztl Z 1991(46): 443-454
028. Gutknecht N, Franzen R, Schippers M, Lampert F: Bactericidal effect of a 980-nm diode- laser in the root canal wall dentin bovine teeth J Clin Laser Med Surg 2004(22): 9-13
029. Kreisler M, Kohnen W, Beck M, Al Haj H, Christoffers AB, Gotz H, Duschner H, Jansen B, D'Hoedt B. Efficacy of NaOCl/H₂O₂ irrigation and GaAlAs laser in decontamination of root canals in vitro. Laser Surg Med 2003; 32(3): 189-196
030. Capodiferro S, Maiorano E, Loiudice AM, Scarpelli F, Favia G: Oral laser surgical pathology: a preliminary study on the clinical advantages of diode laser and on the histopathological features of specimens evaluated by conventional and confocal laser scanning microscopy. Minerva Stomatol 2008(1-2),57: 1-7
031. D'Arcangelo C, Di Nardo Di Maio F, Prosperi GD, Conte E, Baldi M, Caputi S: A preliminary study of healing of diode laser versus scalpel incisions in rat oral tissue: a comparison of clinical, histological, and immunohistochemical results. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2007 Jun;103(6):764-73. Epub 2006 Nov
032. Taylor DL, Schafer SA, Nordquist R, Payton ME, Dickey DT, Bartels KE: Comparison of a high power diode laser with the Nd:YAG laser using in situ wound strength analysis of healing cutaneous incisions. Lasers Surg Med 1997,21(3): 248-254
33. Grigorian AS, Grigoriants LA, Kasparov AS: Experimental- morphological study of the effects of diode laser scalpel with different irradiation parameters (wave length 0.97 mm) on oral mucosa. Stomatologija(Mosk) 2006,85(1): 8-13

034. Stübinger S, Saldamli B, Jürgens P, Ghazal G, Zeilhofer HF: Soft tissue surgery with the diode laser- theoretical and clinical aspects.
Schweiz Monatsschr Zahnmed 2006, 116(8): 812-820
035. Geldi C, Bozkulak O, Tabakoglu HO, Isci S, Kurt A, Gulsoy M: Development of a surgical diode- laser system: controlling the mode of operation.
Photomed Laser Surg 2006,Dec.,24(6): 723-729
036. Gulsoy M, Dereli Z, Tabakoglu HO, Bozkulak O: Closure of skin incisions by 980 nm diode laser welding. Lasers Med Sci 2006, Apr,21(1): 5-10.Epub 2006 Mar 17
037. Saleh HM, Saafan AM: Excision biopsy of tongue lesions by diode laser.
Photomed Laser Surg 2007 Feb, 25(1): 45-49
038. Magid KS, Strauss RA: Laser use for esthetic soft tissue modification.
Dent Clin North Am 2007, Apr, 51(2): 525-545
039. Sano F: Present status and future prospects of laser surgery—advantages of the laser scalpel. Nippon Rinsho 1987, Apr, 45(4): 710-715. PMID : 3302355
040. Barr RE : Laser periodontal treatment and surgical procedures : part1
Pract Proced Aesthet Dent 2004, Nov-Dec,16(10): 747-748. PMID: 15739915
041. Sarver DM: Use of the 810 nm diode laser: soft tissue management and orthodontic applications of innovative technology.
Pract Proced Aesthet Dent 2006, Oct, 18(9): suppl 7-13. PMID: 17139948
042. Walinski CJ: Irritation fibroma removal: a comparison of two laser wavelengths.
Gen Dent 2004, May-Jun, 52(3): 236-238. PMID: 15206255
043. Tamarit- Borrás M, Delgado- Molina E, Berini- Aytes L, Gay- Escoda C:
Removal of hyperplastic lesions of the oral cavity. A retrospective study of 128 cases. Med Oral Patol Oral Cir Bucal 2005, Mar-Apr, 10(2): 151-162.
044. Kravitz ND, Kusnoto B: Soft- tissue lasers in orthodontics: an overview.
Am J Orthod Dentofacial Orthop 2008, Apr, 133(4 Suppl): 110-114
045. Janda P, Sroka R, Mundweil B, Betz CS, Baumgartner R, Leunig A: Comparison of thermal tissue effects induced by contact application of fiber guided laser systems. Lasers Surg Med 2003, 33(2): 93-101

046. Fornaini C, Rocca JP, Bertrand MF, Merigo E, Nammour S, Vescovi P: Nd:YAG and diode laser in the surgical management of soft tissue related to orthodontic treatment. *Photomed Laser Surg* 2007, Oct, 25(5): 381-392
047. Deppe H, Horch HH: Laser applications in oral surgery and implant dentistry. *Lasers Med Sci* 2007, Nov, 22(4): 217-221. Epub 2007 Feb 1
048. Lee EA: Laser assisted gingival tissue procedures in esthetic dentistry. *Pract Proced Aesthet Dent* 2006, Oct, 18(9): suppl 2-6
049. Yeh S, Jain K, Andreana S: Using a diode laser to uncover dental implants in second-stage surgery. *Gen Dent* 2005, Nov- Dec, 53(6): 414-417
050. Romanos GE, Henze M, Banihashemi S, Parsanejad HR, Winckler J, Nentwig GH: Removal of epithelium in periodontal pockets following diode (980 nm) laser application in the animal model: an in vitro study. *Photomed Laser Surg* 2004, Jun, 22(3): 177-183
051. Judy MM, Matthews JL, Aronoff BL, Hults DF: Soft tissue studies with 805 nm diode laser radiation: thermal effects with contact tips and comparison with effects of 1064 nm Nd:YAG laser radiation. *Lasers Surg Med* 1993, 13(5): 528-536
052. Niemz MH: *Laser- Tissue Interactions*. Springer Verlag Berlin, Heidelberg, N.Y. 1996
053. Berlien HP, Müller G: *Angewandte Lasermedizin*. Bd 1. Ecomed Verlagsges. Landsberg , 1999
054. Berlien HP, Müller G: *Angewandte Lasermedizin*. Bd 2. Ecomed Verlagsges. Landsberg , 1999
055. Müller GJ, Ertl T: *Angewandte Laser- Zahnheilkunde* Ecomed Verlagsges. Landsberg/ Lech 1995
056. Vahl J, Van Benthem H: *Laser in der Zahnmedizin*. Quintessenz Verlag Berlin 1992
057. Gutknecht N: *Lasertherapie in der zahnärztlichen Praxis*. Quintessenz Verlag Berlin, Chicago, London 1999.
058. Miserendino LJ, Pick RM: *Lasers in Dentistry*. Quintessenz Verlag Berlin, Chicago 1995.
059. Widgor HA, Walsh Jr JT, Featherstone JDB, Visuri SR, Fried D, Waldvogel JL: *Lasers in Dentistry*. *Laser Surg Med* Vol 16, 103-133, 1995.

060. Hopp M, Schlar N: Laser in der täglichen Anwendung. Teil 2: Fasergestützte Laser Zahn Prax 6: 90-104, 2003
061. Hopp M, Schlar N, Perez-Canto A, Biffar R: Entfernung dystopischen Zungengewebes mittels Diodenlaser (980 nm) aus der Retromolarregion. Z Laserzahnheilkunde 3, 185-189, 2004
062. Hopp M, Bogusch G, Biffar R: Speicheldrüsenexstirpation aus der Unterlippe mittels Diodenlaser (980 nm). Z Laserzahnheilkunde 2, 99-104, 2004
063. Knappe V: Diode Lasers in Berlien H P, Müller G: Applied Laser Medicine. 1. Auflage, Springer Verlag, Berlin, New York, 61-71, 2003
064. Adams TC, Pang PK: Lasers in aesthetic dentistry. Dent Clin North Am 48: 833-860, 2004
065. Coluzzi DJ: Fundamentals of dental lasers: science and instruments. Dent Clin North Am, 48,: 751-770, 2004
066. Jesse K: Laser- Grundlagen und moderne Trends. 1. Auflage VDE Verlag, Berlin- Offenbach, pp 35-54, 1999
067. Meister J, Franzen R, Apel C: Grundlagen der Laserzahnheilkunde, Teil 1: Das Licht. Z Laserzahnheilkunde 1, 57-61, 2004a
068. Meister J, Franzen R, Apel C: Grundlagen der Laserzahnheilkunde, Teil 3: Die Licht- Gewebe-Wechselwirkung. Z Laserzahnheilkunde 3, 199-204, 2004b
069. Niemz MH: Laser- Tissue interactions- Fundamentals and applications. 3.Auflage, Springer Verlag, Berlin- New York, pp 45-149, 2003
070. Walsh LJ: The current status of laser applications in dentistry. Aust Dent J, 48, 146-155, 2003
071. Welch AJ: The thermal response of laser irradiated tissue. IEEE J Quant Electron 20, 1471-1481, 1984
072. Gutknecht N, et al.: Proceedings of the 1st International Workshop of Evidence Based Dentistry on Lasers in Dentistry. Quintessence Publishing Co.Ltd., 2007
073. Duck FA: Physical properties of tissue. Academic Press London- Toronto, 1990
074. Ishii J, Fujita K, Komori T: Laser surgery as a treatment for oral leukoplakia. Oral Oncol 39(8): 759-769, 2003

075. Katzir A: Lasers and optical fibers in medicine.
Academic Press, Inc, San Diego- Boston, 1993
076. Nathan MI, Dumke WP, Burns G, Dill Jr FH, Lasher G: Stimulated emission of radiation from GaAs p-n junctions. Appl Phys Lett 1(3): 62-64, 1962
077. Prahl SA: Optical absorption of hemoglobin.
Oregon Medical Laser Center: <http://omlc.ogi.edu/spectra/hemoglobin/index.html>
078. Strassl M, Kassenbacher A, Wintner E: Ultrashort laser pulses in dentistry.
J Oral Laser Appl 2(4): 213-222, 2002
079. Stabholz A, Zeltzer R, Sela M, Peretz B, Moshonov J, Ziskind D: The use of lasers in dentistry: principles of operation and clinical applications.
Compend 24: 811-824, 2003
080. Kreisler M, Daublander M, Willershausen-Zonnchen B, d`Hoedt B: Effect of diode laser irradiation on the survival rate of gingival fibroblast cell cultures.
Laser Surg Med 28(5): 445-450, 2001a
081. Rastegar S, Jacques SL, Motamedi M, Kim BM: Theoretical analysis of equivalency of high- power diode laser and Nd:YAG laser for coagulation.
Lasers in Dentistry, SPIE 1646: 150-160, 1992
082. Theodoro LH, Haypek P, Bachmann L, Garcia VG, Sampaio JE, Zezell DM, Eduardo C de P: Effect of Er:YAG and diode laser irradiation on the root surface: morphological and thermal analysis. J Periodontol 74(6): 838-843, 2003
083. White JM, Gekelman D, Shin K, Park X, Swenson T, Rouse B, Tran K, Bullard S, Scott- Beckles B, Oto MJSB, Yamamoto A: Laser interaction with dental soft tissues: What do we know from our years of applied Scientific research?
Lasers in Dentistry, SPIE 4610: 39-48, 2002
084. Adams TC, Pang PK: Lasers in Aesthetic Dentistry.
Dent Clin North Am 48: 833-860, 2004
085. Bach G, Neckel C, Mall C, Krekeler G: Conventional versus laser- assisted therapy of periimplantitis: a five year comparative study.
Implant Dent 9: 247-251, 2000
086. Clayman L, Kuo PC: Soft tissue excision techniques. In Clayman L, Kuo P :
Lasers in Maxillofacial Surgery and Dentistry. Thieme Publ., Stuttgart:63-83, 1997

087. Convissar RA, Diamond LB, Fazekas CD: Laser treatment of orthodontically induced gingival hyperplasia. *Gen Dent* 44(1): 47-51, 1996
088. Goharkay K, Moritz A, Wilder –Smith P, Schoop U, Kluger W, Jakolitsch S, Sperr W: Effects on oral soft tissue produced by a diode laser in vitro. *J Oral Surg* 36: 932-937, 1978
089. Olbricht SM, Stern RS, Tang SV, Noe JM, Arndt KA : Complications of cutaneous laser surgery. A survey. *Arch Dermatol* 123(3): 345-349, 1987
090. Pinheiro AL, Limeira Jr FA, Gerbi MEM, Ramalho LMP: Effect of 830 nm Laser Light on the Repair of Bone Defects Grafted with Inorganic Bovine Bone and Decalcified Cortical Osseous Membrane. *J Clin Laser Med Surg* 21(6): 383-388, 2003
091. Romanos GE: Atlas der chirurgischen Laserzahnheilkunde. Urban und Fischer Verlag , München, 1999
092. Romanos GE, Nentwig GH: Present and future of lasers in oral soft tissue surgery: clinical applications. *J Clin Laser Med Surg* 14(4): 179-184, 1996
093. Strauss RA: Laser management of discrete lesions. In : Laser applications in Oral Maxillofacial Surgery. Eds. Catone G A, Alling C C, Saunders, Philadelphia, 1997
094. White JM, Goodis HE, Horton J, Gold S: Current status of lasers in soft tissue dental surgery. *J Periodontol* 65(7): 733-735, 1994
095. Yousuf A, Hossain M, Nakamura Y, Yamada Y, Kinoshita J, Matsumoto K: I Remova of gingival melanin pigmentation with the semiconductor diode laser :. a case report. *J Clin Laser Med Surg* 18: 263-266, 2000
096. Zeinoun T, Nammour S, Dourov N, Aftimos G, Luomanen M: Myofibroblasts in healing laser excision wounds. *Lasers Surg Med* 28: 74-79, 2001
097. Al-Watban FA, Zhang XY: The comparison of effects between pulsed and CW lasers on wound healing. *J Clin Laser Med Surg* 22(1): 15-18, 2004
098. Gutknecht N, Gogswaardt D, Conrads G, Apel C, Schubert C, Lampert F: Diode laser radiation and its bactericidal effect in root canal wall dentin. *J Clin Las Med Surg* 57: 57-60, 2000

099. Parkins F: Lasers in pediatric and adolescent dentistry.
Dent Clin N Amer 44: 821-830, 2000
100. Santaella MR, Braun A, Matson E, Frentzen M: Effect of diode laser and fluoride varnish on initial surface demineralisation of primary dentition enamel : an in vitro study. Int J Paed Dent 14: 199-203, 2004
101. Tuner J, Hode L: The laser therapy handbook. Grängesberg, 231, 2004
102. Gherlone EF, Maiorana C, Grassi RF, Ciancaglini R, Cattoni F: The use of 980 nm diode and 1064 nm Nd:YAG laser for gingival retraction in fixed prostheses. J Oral Laser Applications 4: 183-190, 2004
103. Glauche CE, de Freitas PM, Vieira ND, Marques JL: Qualitative microanalysis of ions and ultrastructural changes in dentin exposed to laser irradiation. Lasers Surg Med 36(4): 334-339, 2005
104. Gutknecht N, Zimmermann R, Lampert F: Lasers in periodontology: state of the art. J Oral Laser Applications 1: 169-179, 2001a
105. Kreisler M, Haj HA, d`Hoedt B: Clinical efficacy of semiconductor laser application as an adjunct to conventional scaling and root planning. Lasers Surg Med 37: 350-355, 2005
106. Mester E, Spiry T, Szende B, Tota JG: Effect of laser rays on wound healing. Am J Surg 22: 532-535, 1971
107. Patruta S, Franz A, Beer F, Goharkhay K, Sperr W, Moritz A: In vitro stimulatory effect of diode laser on the secretion of tissue inhibitor of matrix metalloproteinases 1 (TIMP-1) in human gingival fibroblast: A preliminary report. J Oral Laser Applications 4(3): 167-170, 2004
108. Silveira LB, Silveira JB, Ribeiro MS, Garrocho AA, Novelli MD, Marigo H, Growth EB, et al.: In vivo study on mast cells behavior following low intensity visible and near infrared laser radiation (abstract 304).Lasers Surg Med, 82, 2002
109. Bach G, Koch HK, Hellerich U, Venzke T: Konventionelle Diodenlaser versus Hochpulstechnik. Laser Journal 1/2008: 18-20
110. Deppe H, Horch HH: Laseranwendungen in der chirurgischen Zahn-,Mund- und Kieferheilkunde. Laser Zahnheilkunde 2007(4): 219-222, 2007

111. Horch HH, Deppe H: Laser anwendungen im Mund-Kiefer-Gesichtsbereich. In Horch HH: Mund-Kiefer-Gesichtschirurgie. München: Elsevier, 2007: 797-820
112. Roodenburg JL, Witjes MJ, deVeld DC; Tan IB, Nauta JM: Lasers in dentistry.8. Use of lasers in oral and maxillofacial surgery. Ned Tijdschr Tandheelkd 2002: 109:470
113. Maiorana C, Salina S: Oral soft tissue pathologies: long term evaluation after laser treatment. J Oral Laser Appl 2001(1):20
114. Crippa R: Report on 103 cases of precancerous and other oral lesions treated with a diode laser. J Oral Laser Appl 2001(1): 26
115. Bach G: Er:YAG versus Diode- Zwei Wellenlängen"made in Baden-Württemberg". Laser Zahnheilkunde 2007(4): 283-287
116. Andersen K: Laser technology- a surgical tool of the past, present and future. Aorn J 78: 794-802, 2003
117. Eichler J, Schneeweiss- Wolter C: Wirkung von Laserstrahlung auf Gewebe. MTA 14: 255, 1999
118. McKenzie AL: Physics of thermal processes in laser- tissue interaction. Phys Med Biol 35: 1175-1209, 1990
119. Sullins KE: Diode laser and endoscopic laser surgery. Vet Clin North Am Small Anim Pract. 32: 639-648, 2002
120. Svaasand LO, Boerslid T, Oeveraasen M: Thermal and optical properties of living tissue: application to laser-induced hyperthermia. Lasers Surg Med 5: 589-602, 1985
121. Welch AJ, Motamedi M, Rastegar S, LeCarpentier GL, Jansen D: Laser thermal ablation. Photochem. Photobiol. 53: 815-823, 1991
122. Lippert BM, et al. : Wound healing after laser treatment of oral and oropharyngeal cancer. Lasers Med Sci. 18(1): 36-42, 2003
123. Grossenbacher R: Laser surgery in the oral cavity. J Clin Laser Med Surg 10(1): 19-21, 1992
124. Guttenberg SA, Emery RW: Laser resurfacing and oral and maxillofacial surgeons. Dermatol. Surg. 25(6): 519, 1999

125. Coletton S: Laser treatment of atypical oral pigmentation. A case report.
N Y State Dent J. 72(2): 22-23, 2006
126. Damante CA, et al : Histomorphometric study of the healing of human oral mucosa after gingivoplasty and low-level laser therapy. Lasers Surg Med.35(5): 377-384, 2004
127. Eichler J, Eichler HJ: Laser. Springer Verlag, Berlin- Heidelberg- New York, 1991
128. Manni J: Surgical diode lasers. J Clin Laser Med Surg 10: 377-380, 1992
129. Pick RM, Powel GL: Lasers in dentistry, soft tissue procedures.
Dent Clin North Am 37: 281-297, 1993
130. Scherer H, Fuhre A, Hopf J: Current status of laser surgery of benign diseases in the area of soft palate. Laryngol Rhinol Otol 73: 14-18, 1994
131. Morita T, Okamura Y, Ookub K, et al. : A new instrument for hemostatic vascular management in endoscopic surgery: diode laser(STATLaser- SDL) and it's hand piece. (Dual Hook). J Clin Laser Med Surg 17: 57-61, 1999
132. Claus PN: Comparative study on CW- mode versus pulsed mode in AlGaAs Diode lasers. SPIE Proc 4249: 44-49, 2001
133. Krause F, Frentzen M: Dioden- Laser in der zahnmedizinischen Anwendung. DGL Zeitschrift für Laserzahnheilkunde 3/07: 151-158, 2007
134. Coluzzi DJ: Lasers and soft tissue curettage. An update.
Compend Contin Educ Dent 23: 1104-1111, 2002
135. Dederich DN: Laser curettage: an overview.
Compend Contin Educ Dent 23: 1097-1103, 2002
136. Schwarz F, Sculean A, Berakdar M, Szathmari L, Geoerg T, Becker J: In vivo and in vitro effects of an Er:YAG laser, a diode laser, and scaling and root-planning on periodontally diseased root surfaces: a comparative histologic study. Lasers Surg Med 32: 359-366, 2003
137. Fontana CR, Kurachi C, Mendonca CR, Bagnato VS: Temperature variation at soft periodontal and rat bone tissue during a medium- power diode laser. exposure. Photomed Lasres Surg 22: 519-522, 2004

138. Gold SJ, Vilardi MA: Pulsed laser beam effect on gingiva.
J Clin Periodontol 21(6): 391-396, 1994
139. Borrajo JL, Varela LG, Castro LG, Rodriguez- Nunez I, Torreira MG: Diode laser (980 nm) as an adjunct to scaling and root planning.
Photomed Laser Surg 22: 509-512, 2004
140. Kreisler M, Al Haj H, Daublander M, Gotz H, Duschner H, Willershausen B, D`Hoedt B: Effect of diode laser irradiation on root surface in vitro.
J Clin Laser Med Surg 20: 63-69, 2002
141. Meister J, Franzen R: Dentale Lasersysteme Teil II: Der Dioden- Laser.
DGL Zeitschrift für Laserzahnheilkunde 3/07: 185-189, 2007
142. Alam T, Dawasaz AA, Thukral N, Jangam D: Surgical Diode Laser Excision for Peripheral Cemento-Ossifying Fibroma: A Case Report and Literature Review
J Oral Laser Applications 8: 43-49, 2008
143. Convissar RA: An overview of Laser wavelength used in dentistry.
Dent Clin North Am 44: 753-765, 2000
144. Kreisler M, Meier C, Stender E, Daubländer M, Willershausen- Zönnchen B, de Hoedt B : Effect of diode laser irradiation on the attachment rate of periodontal ligament cells : an in vitro study. J Periodontol Oct, 72(10): 1312-1317, 2001
145. Black P: Der Dentallaser in der oralen Chirurgie-Masterthese-Teil2
Laser Journal 1/2007: 16-21, 2007
146. Van Hillegersberg R: Fundamentals of Laser Surgery
Eur J Surg Jan 1997, 163(1): 3-12, 1997
147. Haina D, Landthaler M, Braun- Falco O, Waidelich W: Comparison of the maximum coagulation depth in human skin for different types of medical lasers.
Lasers Surg Med 7(4): 355-362, 1987
148. Sullins KE: Diode Laser and endoscopic laser surgery.
Vet Clin North Am Small Anim Pract May 32(3): 639-648, 2002
149. McKenzie AL: Physics of thermal processes in laser- tissue interaction.
Phys Med Biol Sep 35(9): 1175-1209, 1990
150. Reinisch L: Laser physics and tissue interaction.
Otolaryngol Clin North Am Dec 29(6): 893-914, 1996

151. Coluzzi DJ: Fundamentals of dental lasers: science and instruments.
Dent Clin North Am 48: 751-770, 2004
152. Niemz MH: Laser- tissue interactions- Fundamentals and applications.
2. Auflage Springer Verlag, Berlin- New York: 45-149, 2003
153. Maiorana C, Salina S: Versatility of Superpulsed Diode Laser in Oral Surgery:
A clinical report. J Oral Laser Appl 6: 193-199, 2006
154. Gold SI, Villardi MA: Pulsed laser beam effects on gingiva.
J Clin Periodontol 21: 391-396, 1994
155. Kamp M: Halbleiterlaser. Vorlesung WS 2006, Technische Physik,
Uni Würzburg, 2006
156. Schäfer O: Die Entwicklung des Diodenlasers bis zur Digitalpulstechnik.
Laser Journal 2/2008: 6-7, 2008
157. Nanami T, Shiba H, Ikeuchi S, Nagai T, Asanami S, Shibata T: Clinical applications and basic studies of laser in dentistry and oral surgery.
Keio J Med Dec.1993, 42(4): 199-201; 1993