

odobleach[®]

Laser power bleaching

elexxion odobleach[®]

Maximum effectiveness for all diode lasers



Image subject to model

elexxion
dental laser

Scope of delivery: everything you need!

This professional treatment set includes:

- 1 **odobleach**® double syringe + 3 mixing tips
- 1 gingival protector + 3 application tips
- 1 after-bleaching care

Each **odobleach**® pack contains everything you need to restore your patients' teeth to gleaming white in just one session. This helps you to plan your material and time and to invoice them accordingly.

Technical data:

- **odobleach**® + elexxion diode laser = maximum effectiveness
- Individually doseable quantities
- Controlled effect
- Whitening by up to nine stages depending on colour type/discolouration
- Efficient application in 30–40 minutes
- Ideal pH value: 6.5

The following dental lasers are already prepared for **odobleach**®:



elexxion **pico**



elexxion **nano**



elexxion **claros**

elexxion **power bleaching**

When it comes to bleaching, two things take centre stage: happy patients and a highly and consistently cost-effective method of treatment.

With the high-quality diode lasers from elexxion, you have access to a range of equipment with functions and a cost/benefit ratio that make it ideal for any practice, from the brand-new, mobile **elexxion pico** with state-of-the-art lithium-ion batteries to the **elexxion claros**®.

Safer and more economical!

Conventional H₂O₂ bleaching gels require long exposure times, often necessitate multiple sessions and produce an acidic pH value with the risk of loss of sensation.

odobleach® is different! With **odobleach**®, titanium dioxide and hydrogen peroxide work in tandem, boosted by the intensity of the elexxion diode laser. **odobleach**® is specially matched to the wavelength of the elexxion equipment.

The result: short exposure times, no damage to the tooth surface, almost no hypersensitive delayed reactions and a high degree of cost-effectiveness.

Interested?

Find out more about the complete elexxion range and **odobleach**® at www.elexxion.com

